

## El INVIMA informa a los usuarios en general que el Grupo de Tecnovigilancia ha emitido una comunicación relacionada con una Alerta asociada a:

<b>NOMBRE DEL DISPOSITIVO MÉDICO</b>	PET/CT PHILIPS - Equipos Combinados de Tomografía Por Emisión de Positrones y Tomografía Computarizada de Rayos X
<b>NO. IDENTIFICACIÓN RISARH</b>	A1402-66
<b>REFERENCIAS DEL DISPOSITIVO MEDICO</b>	GEMINI TF 16 PET/CT y GEMINI TF 64 PET/CT
<b>REGISTRO SANITARIO</b>	INVIMA 2008EBC-0001538
<b>INDICACIONES Y USO ESTABLECIDOS</b>	La tecnología de tomografía por emisión de positrones y la tomografía computarizada de rayos x.
<b>NOMBRE DEL FABRICANTE</b>	PHILIPS Medical Systems (Cleveland) Inc.
<b>DESCRIPCION DEL PROBLEMA</b>	Se informa que el servidor de reconstrucción PET (PRS) puede bloquearse después de la adquisición de CT en bajas dosis y puede no permitir la adquisición de PET lo que resultaría en estudios incompletos, por ende potenciales diagnósticos erróneos y eventos adversos sobre el paciente.
<b>FUENTE</b>	Anexo
<b>FECHA DE NOTIFICACION</b>	14 de Febrero de 2014

### RECOMENDACIÓN:

En caso de identificar la existencia del producto mencionado anteriormente comuníquese con su proveedor quien determinara las acciones que se llevaran a cabo.

Es importante mantener un estado de alerta, realizando un seguimiento permanente a los productos que se fabrican y/o comercializan en el país, divulgando la información de seguridad respectiva entre los profesionales de la salud que realizan uso de estos recursos tecnológicos.

Para mayor información comuníquese al teléfono 2948700 extensión 3880 en Bogotá, ó al correo electrónico [tecnovigilancia@invima.gov.co](mailto:tecnovigilancia@invima.gov.co)

ANEXO

www.ecri.org • Printed from Health Devices Alerts on Friday, February 14, 2014 Page 1

**A21544 01 - High Priority Medical Device Alert**

**Medical Device**  
**Ongoing Action**

Updated: February 11, 2014

**UMDNS Terms:**

- Scanning Systems, Computed Tomography/Positron Emission Tomography [20161]
- Scanning Systems, Positron Emission Tomography [16375]

**Suggested Distribution:**

- Clinical/Biomedical Engineering
- Diagnostic Imaging
- Information Technology

**Geographic Regions:**

- (Impact in specific regions has not been identified or ruled out at the time of this posting)
- Worldwide

**\*Philips—Gemini Positron Emission Tomography/Computed Tomography Systems: PRS Database May Lock Up, Potentially Resulting in Additional Radiation Exposure [Update]**

**Product Identifier:**

[February 7, 2014] Section 1

	Model Nos.:	Serial Nos.:
16-Slice	882470	7003, 7004, 7005, 7006, 7008, 7009, 7010, 7011, 7014, 7015, 7017, 7019, 7020, 7021, 7022, 7024, 7025, 7026, 7031, 7034, 7035, 7036, 7039, 7043, 7045, 7049, 7050, 7052, 7053, 7054, 7058, 7060, 7062, 7064, 7066, 7068, 7069, 7074, 7075, 7077, 7080, 7081, 7084, 7085, 7086, 7087, 7089, 7093, 7094, 7095, 7096, 7097, 7098, 7099, 7108, 7112, 7114, 7115, 7118, 7120, 7121, 7122, 7123, 7124, 7125, 7128, 7129, 7131, 7135, 7137, 7138, 7139, 7140, 7142, 7143, 7145, 7147, 7148, 7156, 7159, 7161, 7163, 7166, 7168, 7173, 7174, 7176, 7177, 7181, 7183, 7184, 7185, 7186, 7187, 7193, 7194, 7195, 7196, 7197, 7198, 7199, 7200, 7207, 7208, 7210, 7211, 7215, 7216, 7223, 7224, 7226, 7228, 7230, 7233, 7235, 7239, 7240, 7242, 7501, 7503, 7507, 7508, 7511, 7514, 7517, 7518, 7520, 7521, 7523, 7527, 7528, 7529, 7530, 7531, 7532, 7533, 7536, 7538, 7540, 7545, 7546, 7554, 7559, 7564, 7566, 7568, 7570, 7573, 7581, 7583, 7584, 7586, 7590, 7591, 7104M, 7106M, 7170M, 7582M

64-Slice	882471	7007, 7012, 7013, 7016, 7027, 7028, 7029, 7030, 7032, 7037, 7042, 7044, 7048, 7055, 7057, 7059, 7061, 7063, 7065, 7071, 7072, 7076, 7078, 7079, 7082, 7083, 7088, 7091, 7092, 7102, 7103, 7105, 7116, 7117, 7119, 7126, 7132, 7136, 7141, 7144, 7146, 7149, 7152, 7153, 7154, 7158, 7160, 7164, 7167, 7169, 7171, 7172, 7175, 7178, 7179, 7180, 7182, 7188, 7189, 7190, 7192, 7201, 7202, 7203, 7204, 7205, 7206, 7209, 7212, 7213, 7214, 7217, 7218, 7219, 7221, 7227, 7229, 7231, 7234, 7236, 7238, 7241, 7245, 7251, 7252, 7253, 7254, 7255, 7256, 7257, 7258, 7504, 7506, 7509, 7510, 7512, 7515, 7516, 7519, 7524, 7525, 7535, 7537, 7539, 7542, 7543, 7544, 7548, 7549, 7550, 7552, 7553, 7555, 7556, 7557, 7558, 7560, 7561, 7562, 7563, 7567, 7569, 7572, 7574, 7575, 7577, 7578, 7579, 7585, 7587, 7589, 7592
----------	--------	--

Big Bore	882476	9004, 9005, 9006, 9008, 9009, 9010, 9011, 9012, 9013, 9014, 9015, 9016, 9017, 9018, 9019, 9020, 9021, 9022, 9023, 9201, 9202, 9203, 9204, 9205, 9206, 9207, 9208, 9209, 9210, 9211, 9213, 9214, 9215, 9216, 9217, 9218, 9219, 9220, 9221, 9222, 9223, 9224, 9225, 9226, 9228, 9229, 9231
LXL	882412	14701, 14702, 14704, 14705, 14706, 14707, 14708, 14709, 14710, 14711
Ready	882473	7040, 7048, 7070, 7095, 7111, 7113, 7127, 7130, 7134, 7150, 7151, 7155, 7157, 7162, 7237, 7249

TruFlight Select	882438	17002, 17003, 17004, 17005, 17006, 17007, 17008, 17009, 17010, 17011, 17012, 17013, 17014, 17015, 17017, 17018, 17019, 17021, 17022, 17023, 17024, 17025, 17026, 17028
------------------	--------	--

[Capital Equipment]

[January 2, 2014] Section 2

Gemini TF Positron Emission Tomography (PET)/Computed Tomography (CT) Systems:	Software Versions:
(1) 16-Slice	(1 and 2) 3.5.1, 3.5.1.1, 3.5.2, 3.5.2.1
(2) 64-Slice	
(3) Big Bore	(3) 3.6, 3.6.1, 3.6.2
(4) LXL	(4 and 5) 3.5.2
(5) Ready	
(6) TruFlight Select	(6) V3.5.3

[Capital Equipment]

Manufacturer:

- Philips Healthcare North America 3000 Minuteman Rd, Andover, MA 01810-1099, United States

**Summary:** This Alert provides additional information based on FDA's Center for Devices and Radiological Health (CDRH) source material regarding [Alert Accession Number A21544](#). FDA's CDRH states that the manufacturer initiated a recall by Field Safety Notice letter dated December 4, 2013. Additional information is provided in the following field:

- Product Identifier (see Section 1)

**Problem:**

[December 31, 2013]

In a December 4, 2013, Urgent Medical Device Correction letter submitted by an ECRI Institute member hospital and posted by the U.K. Medicines and Healthcare Products Regulatory Agency (MHRA), Philips states that the PET reconstruction server (PRS) database in the above systems may lock up after the low-dose CT acquisition is complete and may not allow the PET acquisition to begin, resulting in an incomplete study. Philips states that if this problem occurs, users cannot fix it and a service interaction is required before the system will become operational. Delay of a patient's study may necessitate a rescans and/or patient reinjection if a short half-life isotope or radiopharmaceutical is used, exposing the patient to additional radiation. The manufacturer has not confirmed the information provided in the source material.

**Action Needed:**

The following actions are those listed in [Alert Accession Number A21544](#). Identify any affected product in your inventory by clicking the help button and "Product Info" tab to display the software version. If you have affected product, verify that you have received the December 4, 2013, Urgent Medical Device Correction letter from Philips. A Philips service engineer will contact your facility to arrange to install a software update on affected systems.

**For Further Information:**

North America  
Philips customer care solutions center  
Tel.: (800) 722-9377 (select option 5)  
Europe  
Philips UK Customer Care Service Centre  
Tel.: (0870) 5329741  
Website: [Click here](#).

**References:**

- Great Britain. Medicines and Healthcare Products Regulatory Agency. Philips Medical Systems. Nuclear medicine, PET-CT PET-CT system. GEMINI TF 16, GEMINI TF 64, GEMINI TF Big Bore, TruFlight Select, GEMINI TF Ready and GEMINI LXL PET/CT systems [online]. London: Department of Health; 2013 Dec 20 [cited 2013 Dec 31]. (Field safety notice; reference no. 2013/012/017/291/009). Available from Internet: [Click here](#).
- United States. Food and Drug Administration. Center for Devices and Radiological Health. Medical device recalls: Class 2 recall—Gemini LXL Gemini TF PET/CT16 Gemini TF PET/CT 64 [online]. 2014 Feb 4 [cited 2014 Feb 7]. Available from Internet: [Click here](#).

**Comment:**

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

**Source(s):**

- 2014 Feb 7. FDA CDRH Database. Class II. Z-0932-2014
- 2014 Feb 7. Member Hospital. FSN 88200465-471
- 2014 Feb 7. MHRA FSN. 2013/012/017/291/009