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## ANEXO 1

www.ecri.org • Printed from *Health Devices Alerts* on Friday, October 03, 2014 Page 1

### H0240 - High Priority Medical Device Alert

#### Medical Device Hazard Report

Updated: October 2, 2014

#### UMDNS Terms:

- Ventilators, Intensive Care, Neonatal/Pediatric [14361]
- Ventilators, Intensive Care [17429]

#### Suggested Distribution:

- Pulmonology/Respiratory Therapy
- Critical Care
- NICU
- Clinical/Biomedical Engineering

#### Geographic Regions:

- Worldwide

#### Covidien—Puritan Bennett 980 Ventilators: Graphical User Interface (GUI) May Become Dim if Capacitor Fails

##### Product Identifier:

Covidien Puritan Bennett 980 Ventilators [Capital Equipment]

##### Manufacturer:

- Covidien (US) 15 Hampshire St, Mansfield, MA 02048, United States

##### Problem:

During operation, the graphical user interface (GUI) display for the Covidien Puritan Bennett 980 ventilator may emit an electrical burning odor and become dim. Covidien states that ventilation will continue at all previously set parameters and that although the display is dim, machine settings, patient parameters, and alarms are still displayed. During such an event, a substitute ventilator should be employed.

##### Suppliers Corrective Action

Covidien reports that the problem is caused by a capacitor failure on the GUI backlight inverter board. If this occurs, the display illumination will decrease by 50%, but a duplicate circuit is unaffected and continues to power 2 GUI backlights. Covidien will be issuing a customer notification discussing this problem and an unrelated software update in the near future. Covidien plans to issue the notification by October 3rd depending on health authority communications. Covidien also indicated that field service staff will update the 980 ventilators with a new backlight inverter board.

##### ECRI Institute Recommendations:

Until the inverter board is installed, ECRI Institute recommends that users of the above ventilators do the following:

- Ensure that a substitute ventilator is available in the event that the GUI backlight inverter board fails.
- Always have a manual resuscitator available for manual ventilation in the event of a ventilator failure.

##### Comment:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

##### Source(s):

- 2014 Oct 1.
- ECRI Institute researched member report.