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**ANEXO 1**

www.ecri.org • Printed from *Health Devices Alerts* on Friday, April 11, 2014 Page 1

**A22126 - High Priority Medical Device Alert**

**Medical Device**

**Ongoing Action**

Updated: April 9, 2014

**UMDNS Terms:**

- Workstations, Picture Archiving and Communication System [21968]

**Suggested Distribution:**

- Clinical/Biomedical Engineering
- Diagnostic Imaging
- Information Technology

**Geographic Regions:**

- Worldwide

**GE—Centricity PACS RA1000 Workstations: Previously Displayed Images May Be Rejected in a Reject Image Workflow**

**Product Identifier:**

Centricity Picture Archiving and Communication Systems (PACS) RA1000 Workstations [*Capital Equipment*]  
Software Versions: 3.2, 3.2.0.1, 3.2.0.2, 3.2.1, 3.2.2, 3.2.2.1, 3.2.2.2, 3.2.2.3, 3.2.6, 3.2 SP7, 3.2 SP8, 4.0, 4.0.1, 4.0 SP3

**Manufacturer:**

- GE Healthcare IT540 W Northwest Hwy, Barrington, IL 60010, United States

**Problem:**

In an April 3, 2014, Urgent Medical Device Correction letter submitted by ECRI Institute member hospitals, GE states that in a reject image workflow, in which multiple studies are simultaneously opened on the above workstations and images are rejected from the active displayed study, images from the previously displayed study may also be rejected. GE states that the problem can occur with a previously displayed study that remains open in the background under the following circumstances:

- If one or more, but not all, images in the exam were inadvertently rejected, the image count in the work modes palette identifies fewer images; however, all the images remain viewable.
- For newly acquired exams, if the entire exam was inadvertently rejected, the workstation will return a message stating "Unable to fetch this exam: No images in the exam." The images would then need to be resent to PACS from the modality.
- For historical exams, if the entire exam was inadvertently rejected, the workstation will return a message stating "Unable to fetch this exam: No images in the exam," and the exam is not available on the modality, a GE Healthcare service representative would need to be contacted to restore the exam on the database.

GE also states that there is a DelRejlmj task in the above systems, that if enabled and running, in conjunction with the inadvertently rejected image problem, poses a potential patient safety hazard because of data loss. If the DelRejlmj task runs, the inadvertently rejected images will be permanently deleted from the system and may result in a missed intervention or misdiagnosis. If you do not have the DelRejlmj task running, there is no potential hazard since the inadvertently rejected images remain on the system. GE states that it has received no reports of patient injuries related to the above problems.

**Action Needed:**

Identify any affected product in your inventory. If you have affected product, verify that you have received the April 3, 2014, Urgent Medical Device Correction letter from GE. A GE Healthcare service representative will contact your facility to arrange to correct affected systems at no cost. Until the software update is installed, GE recommends that you perform the following steps:

- (1) Check to see if the DelRejlmj task is active and running.
  - (a) Log into the CA Tool.
  - (b) From the Navigation menu, select "System Monitor."
  - (c) From the pull down menu, select "Dbase Maintenance."
  - (d) Scroll to the DelRejlmj task (this list is alphabetical).
  - (e) Review the Task Status. Possible statuses are "IDLE" or "SUSPEND."
  - (f) If the task status is set to "SUSPEND," the task is not running.
  - (g) If the task is set to "IDLE," move to Step 2.
- (2) Check the Minimum Rejected Image Deletion Time.
  - (a) From the Navigation menu, select "Sys. Configuration."
  - (b) From the pull down menu, select "Pacs Configuration."
  - (c) Scroll to "Minimum Rejected Image Deletion Time" (this list is alphabetical).
  - (d) Review the Value field.
  - (e) If the value field is a negative integer (such as the default setting of -1), the Delete Rejected Image process is NOT running and no patient safety hazard exists. If the

value field is a positive integer (such as 1 or higher), the Delete Rejected Image process is running and a potential patient safety hazard exists.  
(3) A GE Healthcare service representative will remove the DelRejImg task from the task table in the Centricity PACS database. This will prevent the DelRejImg task from running and prevent any rejected images from being deleted. This feature is being removed entirely from the 3.2 release family of Centricity PACS and will be reinstated in a future 4.X workstation release of Centricity PACS.  
(4) Perform the reject image workflow with only one patient study open at a time.

**For Further Information:**

**Americas**  
Tel.: (800) 437-1171  
**Europe, Middle East, Africa**  
Tel.: 33 (1) 30831300  
**Singapore**  
Tel.: 63880932  
**Malaysia**  
Tel.: 1800883911  
**Australia**  
Tel.: 1 (800) 659465  
**New Zealand**  
Tel.: 800659465  
**India**  
Tel.: 1800114567  
**Korea**  
Tel.: 15446119  
**Japan**  
Tel.: 120055919  
**Taiwan**  
Tel.: 800021770  
**China**  
Tel.: 8008108188  
**Hong Kong**  
Tel.: 852-21006288  
Website: [Click here](#)

**Comment:**

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

**Source(s):**

- 2014 Apr 8. Member Hospital. GE letter submitted by ECRI Institute member hospital; Reference No. 85324
- 2014 Apr 9. Manufacturer. Manufacturer confirmed the information provided in the source material.