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ANEXO 1

www.ecri.org • Printed from *Health Devices Alerts* on Monday, August 25, 2014 Page 1

A22871 - High Priority Medical Device Alert

**Medical Device
Ongoing Action**

Updated: August 14,
2014

UMDNS Terms:

- Scanning Systems, Computed Tomography/Single Photon Computed Tomography [24013]

Suggested Distribution:

- Clinical/Biomedical Engineering
- Diagnostic Imaging
- Nuclear Medicine

Geographic Regions:

- Worldwide

Siemens—Symbia Intevo and Symbia T SPECT/CT Systems: Rubber Isolator May Break and Cause Loss of Support for Rotate Motor

Product Identifier:

Symbia Single Photon Computed Tomography/Computed Tomography (SPECT/CT) Systems: (1) Intevo, (2) T [*Capital Equipment*]

Manufacturer:

- Siemens Medical Solutions USA Inc685 East Middlefield Rd, Mountain View, CA 94043-7393, United States

Problem: In a Customer Safety Advisory Notice letter submitted by ECRI Institute member hospitals, Siemens states that the rubber isolator inside the rotate motor mounting assembly of the above systems may degrade and break, potentially causing loss of support for the rotate motor. Siemens also states that this problem may lead to a rupture in the drive belt, potentially damaging the gantry covers and making contact with the patient. Siemens further states that it has received no reports of injury related to this problem.

Action Needed:

Identify any affected systems in your inventory. If you have affected systems, verify that you have received the Customer Safety Advisory Notice letter from Siemens. Siemens states that you may continue to use the above systems; however, if you observe abnormal noise or vibration, discontinue use of the system and notify your local service representative. Your Siemens local service engineer will contact your facility to arrange for mounting plate replacement. Retain a copy of the letter in the system's instructions for use. Forward a copy of the letter to any facility to which you have further distributed affected product, and notify Siemens of the transfer.

For Further Information:

Siemens
America
Tel.: (800) 888-7436
Europe, Middle East, and Africa
Tel.: 49 (9131) 9404000
Asia and Australia
Tel.: 86 (21) 38112121
Website: [Click here](#)

Comment:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2014 Aug 11. Member Hospital. Siemens letter submitted by an ECRI Institute member hospital.
- 2014 Aug 12. Manufacturer. The manufacturer confirmed the information in the source material.
- 2014 Aug 13. BfArM (Germany). 4594/14
- 2014 Aug 13. BfArM (Germany).