**PROSPERIDAD  
PARA TODOS**

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ANEXO 1

www.ecri.org • Printed from Health Devices Alerts on Thursday, February 13, 2014 Page 1

**A21652 - Normal Priority Medical Device Alert**

**Medical Device**

**Ongoing Action**

**Updated:** February 5, 2014

**UMDNS Terms:**

- Analyzers, Laboratory, Clinical Chemistry [15551]

**Suggested Distribution:**

- Clinical Laboratory/Pathology Engineering
- Information Technology

**Geographic Regions:**

- Argentina
- Australia
- Austria
- Belarus
- Brazil
- Bulgaria
- Canada
- Colombia
- Cyprus
- Czech Republic
- Dominican Republic
- Egypt
- Finland
- France
- Germany
- Greece
- Hungary
- India
- Italy
- Kazakhstan
- Latvia
- Lebanon
- Lithuania
- Norway
- Peru
- Poland
- Portugal
- Romania
- Russia
- Serbia
- Slovakia
- South Africa
- South Korea
- Spain
- Sweden
- Turkey
- U.K.
- United Arab Emirates

**Siemens—ADVIA 1800 and 2400 Chemistry Systems: Electrodes May Fail Prematurely during Calibration**

**Product Identifier:**

ADVIA Chemistry Systems: (1) 1800, (2) 2400 [Capital Equipment, Consumable] Firmware Version 2.23

**Manufacturer:**

- Siemens Healthcare Diagnostics Inc 511 Benedict Ave, Tarrytown, NY 10591, United States

**Problem:**

In a January 2014 Urgent Field Safety Notice letter posted by the U.K. Medicines and Healthcare Products Regulatory Agency (MHRA), Siemens states that the above chemistry systems may fail prematurely during calibration because of a chloride bias parameter error. Siemens also states that the operator may receive an NG flag, potentially requiring immediate action.

**Action Needed:**

Identify any affected product in your inventory. If you have affected product, verify that you have received the January 2014 Urgent Field Safety Notice letter and field correction effectiveness check form from Siemens. Your Siemens representative will contact you to arrange for firmware to be reverted back to version 2.22. Review the Urgent Field Safety Notice letter with your facility's laboratory medical director. Complete the field correction effectiveness check form, and return it to Siemens using the information provided on the form. Keep a copy of the Urgent Field Safety Notice letter with your laboratory records, and forward a copy of the letter to any facilities to which you have further distributed affected product.

**For Further Information:**

Siemens technical solutions center

Website: [Click here](#)

**References:**

- Great Britain. Medicines and Healthcare Products Regulatory Agency. Siemens Healthcare Diagnostics. IVDs, clinical chemistry instrumentation. Chloride electrode/ISE firmware V2.23 [online]. London: Department of Health; 2014 Feb 3 [cited 2014 Feb 4]. (Field safety notice; reference no. 2014/001/023/601/011). Available from Internet: [Click here](#).

**Comment:**

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

**Source(s):**

- 2014 Feb 4. MHRA FSN. 2014/001/023/601/011
- 2014 Feb 4. MHRA FSN.
- 10817232 (includes reply form)
- 2014 Feb 5. Manufacturer. The manufacturer confirmed the information provided in the source material.