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ANEXO

www.ecri.org • Printed from Health Devices Alerts on Friday, February 14, 2014 Page 1

A21770 - High Priority Medical Device Alert

**Medical Device
Ongoing Action**

Updated: February 13, 2014

UMDNS Terms:

- Stereotactic Systems, Image-Guided, Radiosurgical, Linear Accelerator [26629]

Suggested Distribution:

- Clinical/Biomedical Engineering
- Radiation Oncology/Medical Physics
- Facilities/Building Management

Geographic Regions:

- Worldwide

Accuray—Synchrony Boom Arm Mounting Assembly Used with CyberKnife Systems: Assembly Joint May Detach and Fall

Product Identifier:

Synchrony Arm Mounting Assembly used with CyberKnife Systems [Capital Equipment]

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302 units distributed

Manufacturer:

- Accuray Inc1240 Deming Way, Madison, WI, 53717, United States

Problem:

In a January 14, 2014, Urgent Device Correction letter posted by the U.K. Medicines and Healthcare Products Regulatory Agency (MHRA), Accuray states that a joint on the boom arm mounting assembly for the above cameras may detach and cause the camera to fall, potentially resulting in serious injury to the patient or user. FDA's center for Devices and Radiological Health (CDRH) states that the manufacturer initiated a correction by Urgent Medical Device Correction letter dated January 14, 2014.

Action Needed:

Identify any affected product in your facility. If you have affected product, verify that you have received the January 14, 2014, Urgent Device Correction letter and acknowledgment of notification form from Accuray. Complete and return the acknowledgement of notification form using the instructions provided in the form. Inspect the joints of the boom arm mounting assembly using the instructions provided in the [letter](#). If any twisting motion or loose connections are suspected in the threaded joints of the boom arm mounting assembly, call the Accuray customer support department. If the boom arm mounting assembly is secure, no other immediate action

is required. Accuray will contact your facility to arrange to secure the mounting assembly.

For Further Information:

Accuray
U.S.
Tel.: (877) 668-8667 (USA)
Outside the U.S.
Tel.: (408) 716-4700
Website: [Click here](#)

References:

- Great Britain. Medicines and Healthcare Products Regulatory Agency. Accuray Incorporated. Radiotherapy Cyber Knife. CyberKnife robotic radiosurgery system [online]. London: Department of Health; 2014 Jan 31 [cited 2014 Feb 11]. (Field safety notice; reference no. 2014/001/027/081/020). Available from Internet: [Click here](#).
- United States. Food and Drug Administration. Center for Devices and Radiological Health. Medical Device Recalls: Class 2 recall—Accuray CyberKnife robotic radiosurgery system [online]. 2014 Jan 27 [cited 2014 Feb 11]. Available from Internet: [Click here](#).

Comment:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2014 Feb 10. MHRA FSN. MHRA Reference No.: 2014/001/027/081/020
- 2014 Feb 10. FDA CDRH Database. Class II. Z-0828-2014
- 2014 Feb 10. MHRA FSN. Accuray letter submitted to MHRA (includes reply form)
- 2014 Feb 12. Manufacturer. The manufacturer confirmed the information provided in the source material.