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**ANEXO 1**

www.ecri.org • Printed from Health Devices Alerts on Friday, July 18, 2014 Page 1

**A22655 - High Priority Medical Device Alert**

**Medical Device  
Ongoing Action**

Updated: July 8, 2014

**UMDNS Terms:**

- Beds, Electric [10347]

**Suggested Distribution:**

- Clinical/Biomedical Engineering
- Critical Care
- Facilities/Building Management
- Home Care
- OR/Surgery

**Geographic Regions:**

- Africa
- Argentina
- Australia
- Brazil
- Canada
- Chile
- China
- Europe
- Japan
- Korea
- Latin America
- Mexico
- Middle East
- Singapore
- South Africa
- Taiwan
- U.K.
- U.S.

**Stryker—InTouch Critical Care Beds: Casters May Break if Brake System Is Activated While Bed Is in Motion**

**Product Identifier:**

InTouch Critical Care Beds [Capital Equipment]  
Model Nos.: 2130, 2131, 2140, 2141, 2151, 2152, 2155, 2156

For affected serial numbers, see the Serial Number Listing sent to your facility.

**Manufacturer:**

- Stryker Medical 3800 E Centre Ave, Portage, MI 49002, United States

**Problem:**

In a June 30, 2014, Urgent Medical Device Notification letter submitted by ECRI Institute member hospitals, Stryker states that it has received reports of incidents in which a caster on the above beds has broken on top of the caster stem where the brake rod inserts into the caster. Stryker states that these broken casters are identified when the beds are hoisted or lifted from the floor during regular maintenance. If more than one caster on a unit is damaged, brake functions may be affected and result in reduced brake holding force. Stryker also states that if only one caster on the bed is damaged, brake functionality continues to meet specifications. Stryker further states that the root cause of this problem is attributable to a potential use error in which the bed's electric brakes are activated while the bed is in motion.

**Action Needed:**

Identify any affected product in your inventory. If you have affected product, verify that you have received the June 30, 2014, Urgent Medical Device Notification letter, operation instructions, brake labels and application instructions, and business reply form from Stryker. Stryker states that the caster braking system on the above beds is designed to be activated when the bed is in a stationary position, as stated in the operations manual. When transporting the bed, the bed should be brought to a complete stop before applying the bed's caster brake system. The caster brake system is intended to keep the bed from moving while in a stationary position. Complete the business reply form, and return it to Stryker using the instructions on the form. Affix the brake labels to the bed control panel using the instructions provided with the labels.

**For Further Information:**

Stryker customer service department  
Tel.: (800) 327-0770, 8 am. to 6 p.m. Eastern time, Monday through Friday  
Website: [Click here](#)

**Comment:**

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

**Source(s):**

- 2014 Jul 7. Member Hospital. Stryker letter submitted by ECRI Institute member hospitals (includes reply form)
- 2014 Jul 8. Manufacturer. Manufacturer confirmed information