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ANEXO 1

www.ecri.org • Printed from Health Devices Alerts on Friday, November 21, 2014 Page 1

A23348 - Normal Priority Medical Device Alert

**Medical Device
Ongoing Action**

Updated: November 20, 2014

UMDNS Terms:

- Viscoelastic Solutions, Intraocular [29717]
- Viscoelastic Solutions [33785]

Suggested Distribution:

- Materials Management
- Ophthalmology

Geographic Regions:

- Worldwide

Abbott Medical Optics (AMO)—Healon EndoCoat Ophthalmic Viscosurgical Devices and Healon Duet Dual Packs: Manufacturer Reminds Users of Proper Syringe Technique to Avoid Finger Grip Dislodgement

Product Identifier:

Healon EndoCoat Ophthalmic Viscosurgical Devices (OVDs) [Consumable]
U.S. Product Nos.: VT585U, 57502000 (packaged in Healon Duet Dual Pack Product No. 10290080)
International Product Nos.: VT585, 57550500 (packaged in Healon Duet Dual Pack Product Nos. 10220010, 10220011, 10220012)

Manufacturer:

- Abbott Medical Optics 1700 E. St. Andrew Place, Santa Ana, CA 92705, United States

Problem: In a November 4, 2014, Urgent Physician Advisory Notice letter submitted by ECRI Institute member hospitals, Abbott Medical Optics (AMO) states that it has received reports of incidents in which the finger grip on the above ophthalmic viscosurgical devices (OVDs) dislodged while the user depressed the syringe plunger. If the finger grip dislodges, the surgeon's hand could accidentally move the cannula further than intended, potentially resulting in trauma to the patient's eye. AMO states that the finger grip may become detached if the user fails to follow the directions for use (DFU) by using a cannula not provided with the device or by using excessive force.

Action Needed:

Identify affected product in your inventory. If you have affected product, verify that you have received the November 4, 2014, Urgent Physician Advisory Notice letter and additional diagram from AMO. AMO will add the new diagram to the affected product's DFU. AMO also recommends the following:

- Only use the 25 gauge cannula provided with affected product.
- Avoid excessive force when using the syringe plunger.
- Check proper function by holding the syringe barrel and gently depressing the plunger rod until the OVD appears at the cannula tip.

Inform all relevant personnel of the information in the Urgent Physician Advisory Notice letter, and provide a copy of the letter to any facility to which you have distributed affected product.

For Further Information:

AMO
(877) 266-4543
Website: [Click here](#)

Comment:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2014 Nov 13. Member Hospital. Abbott letter submitted by ECRI Institute member hospital
- 2014 Nov 20. Manufacturer. Manufacturer confirmed information.