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ANEXO 1

www.ecri.org • Printed from Health Devices Alerts on Monday, September 08, 2014 Page 1

A22982 - High Priority Medical Device Alert

**Medical Device
Ongoing Action**

Updated: September 3,
2014

UMDNS Terms:

- Scanning Systems,
Gamma Camera,
Single Photon
Emission Tomography
[18444]

Suggested Distribution:

- Clinical/Biomedical
Engineering
- Diagnostic Imaging

Geographic Regions:

- Worldwide

**Siemens—e.cam Single-Photon Emission Computed
Tomography Patient Handling Systems: Patients May Be
Injured during System Movement**

Product Identifier:

e.cam Single-Photon Emission Computed Tomography (SPECT) Patient Handling
Systems (PHS) [Capital Equipment]

Manufacturer:

- Siemens Medical Solutions USA Inc2501 N Barrington Rd, Hoffman Estates, IL
60192, United States

Problem: In an April 8, 2014, Customer Safety Advisory Notice letter submitted by an ECRI Institute member hospital, Siemens states that it has received one report of an incident in which a patient was reportedly grasping the integrated armrest of the above systems as the technician pulled the PHS pallet to the home position; during this movement, when the pallet was retracted, the patient's finger was caught between the moving and the stationary parts of the PHS and injured.

Action Needed:

Identify any affected product in your inventory. If you have affected product, verify that you have received the April 8, 2014, Customer Safety Advisory Notice letter from Siemens. Siemens recommends that you monitor the patient at all times while he/she is on the system and during initiation of any system movement, as described in the device operator's manual and Patient Bed Safety Addendum of the operator's manual. Siemens also recommends that you observe the following general good practices:

- Never leave a patient unattended on the system.
- When moving a patient into or out of the gantry, make sure the patient's arms, fingers, hair, and loose items such as clothing are completely positioned on the PHS.
- When moving any components of the system, all areas should be checked to ensure that there are no obstacles in the path of motion.

Retain a copy of the letter in your device operator's manual, and notify all relevant personnel at your facility of the information in the letter. Forward a copy of the letter to any facility to which you have further distributed affected product, and inform Siemens of the new owner.

For Further Information:

Siemens

Tel.: (800) 888-7436, 7 am. to 8 p.m. Eastern time, Monday through Friday

Website: [Click here](#)

Comment:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2014 Sep 3. Member Hospital. Siemens letter submitted by an ECRI Institute member hospital.
- 2014 Sep 3. Manufacturer. The manufacturer confirmed the information in the source material.