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ANEXO 1

www.ecri.org • Printed from Health Devices Alerts on Friday, September 19, 2014 Page 1

A22998 - Normal Priority Medical Device Alert

**Medical Device
Ongoing Action**

Updated: September 16, 2014

UMDNS Terms:

- Analyzers, Laboratory, Clinical Chemistry/Immunoassay [20821]

Suggested Distribution:

- Clinical/Biomedical Engineering
- Clinical Laboratory/Pathology

Geographic Regions:

- Canada
- China
- Colombia
- Croatia
- Czech Republic
- France
- Germany
- India
- Italy
- Japan
- Korea
- Kuwait
- Malaysia
- Mexico
- Philippines
- Puerto Rico
- Russia
- Singapore
- Slovakia
- South Africa
- Spain
- Sweden
- Switzerland
- Taiwan
- Thailand
- Turkey
- U.K.
- U.S.
- Vietnam

Beckman Coulter—UniCel DxI Access Immunoassay Systems and UniCel DxC Synchron Access Clinical Systems: Nonconforming Sample Probe Cables May Have Been Installed

Product Identifier:

UniCel DxI Access Immunoassay Systems: (1) 600, (2) 800, and UniCel DxC Synchron Access Clinical Systems: (3) DxC 660i, (4) DxC 680i, (5) DxC 860i, (6) 880i [Capital Equipment]

Manufacturer:

- Beckman Coulter (UK) Ltd Oakley Court Kingsmead Business Park, High Wycombe HP11 1JU, England

Problem:

In an August 29, 2014, Urgent Field Safety Notice letter posted by the U.K. Medicines and Healthcare Products Regulatory Agency (MHRA), Beckman Coulter states that if the above systems are newly manufactured and/or have recently been serviced, they may contain nonconforming sample probe cables. Beckman Coulter also states that use of a nonconforming sample probe cable may cause a delay in reporting results; however, this will not cause erroneous results. Beckman Coulter further states that a nonconforming sample probe may also cause level sense errors, which may prevent tests from starting or automatically cancel an in-process test.

Action Needed:

Identify any affected instruments in your inventory. Only newly manufactured instruments shipped from Beckman Coulter between May 22 and August 8, 2014, or instruments serviced during the same time period, are potentially affected by this problem. If you have affected instruments, verify that you have received the August 29, 2014, Urgent Field Safety Notice letter and Response Form from Beckman Coulter. Beckman Coulter states that you may continue to operate affected instruments normally and that your Beckman Coulter service representative will contact your facility to schedule a service visit to ensure that a conforming sample probe cable is installed on your system(s). Complete the Response Form, and return it to Beckman Coulter using the instructions on the form. Retain a copy of the Urgent Field Safety Notice letter as part of your laboratory Quality System documentation, notify all relevant personnel at your facility of the information in the letter, and forward a copy of the letter to any facility to which you have further distributed affected product.

For Further Information:

Beckman Coulter local representative or Beckman Coulter customer support hotline
Tel.: (0845) 6001345
E-mail: techsupportuk@beckman.com
Website: [Click here](#)

References:

- Great Britain. Medicines and Healthcare Products Regulatory Agency. Beckman Coulter. IVDs, clinical chemistry assay kit immunoassay. Sample probe cable—for use with the UniCel DxI family of immunoassay systems [online]. London: Department of Health; 2014 Sep 8 [cited 2014 Sep 9]. (Field safety notice; reference no. 2014/008/029/601/003). Available from Internet: [Click here](#).

Comment:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this

alert, refer to our [HDA Format Guide](#).

Source(s):

- 2014 Sep 9. MHRA FSN. 2014/008/029/601/003
- 2014 Sep 9. MHRA FSN. FSN 23515
- 2014 Sep 16. Manufacturer. Beckman Coulter confirmed the information provided in the source material.