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ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Wednesday, April 22, 2015 Page 1

[Normal Priority] - S0273 : Philips—BrightView Gamma Cameras: Detectors May Be Damaged by Hitting the Floor during Repositioning [ECRI Institute User Experience Network] Medical Device Special Report

Published: Monday, March 30, 2015

UMDNS Terms:

- Scanning Systems, Gamma Camera, Single Photon Emission Tomography [18444]
- Scanning Systems, Computed Tomography/Single Photon Computed Tomography [24013]

Product Identifier: Single Photon Emission Computed Tomography (SPECT) Gamma Cameras: (1) Brightview, (2) Brightview X, (3) Brightview XCT [*Capital Equipment*]

Geographic Regions: Worldwide

Manufacturer(s): Philips Healthcare North America 3000 Minuteman Rd, Andover, MA 01810-1099, United States

Suggested Distribution: Clinical/Biomedical Engineering, Diagnostic Imaging, Nuclear Medicine

Problem:

1. The detectors on the above gamma camera systems may come into contact with the floor while being repositioned between scans when the system is in the preprogrammed motion mode.
2. This problem may result in serious damage to the gamma camera.

ECRI Institute Recommendations:

In addition to the manufacturer's instructions and the existing hospital's procedures, ECRI Institute recommends the following:

Nuclear Medicine Technologists/Operators:

1. Ensure that the detectors are raised to the isocenter of the gantry before repositioning the detectors for the next study.
2. Observe the camera during repositioning while staying within reach of the emergency stop switch. Activate the emergency stop switch if contact with the floor or any other objects is likely.
3. Report any incidents of the detector contacting the floor to ECRI Institute and the manufacturer.

Manufacturer Perspectives

Philips has confirmed that they have replicated the problem and have not taken corrective action at this time.

ECRI Institute Perspectives

ECRI Institute believes gamma cameras should be designed to prevent the detector from hitting the floor during automatically controlled motion.

Background:

1. A member hospital reported that while repositioning the detectors on a Philips BrightView gamma camera using the preprogrammed motion (PPM) mode, one of the detectors collided with the floor.
2. The camera had just completed a pediatric standing study, for which the detectors were positioned below the isocenter of the system and close to the floor.
3. The detectors were being repositioned in preparation for the next patient, while there was no patient in the scan room.
4. The PPM mode was used to complete all the repositioning steps.
5. The detector came into contact with the floor and was damaged.

©2015 ECRI Institute
5200 Butler Pike, Plymouth Meeting, PA 19462-1298, USA
May be reproduced by subscribing institution for internal distribution only.