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ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Thursday, August 20, 2015 Page 1

[High Priority] - A24806 : Siemens—RAPIDPoint 500 Blood Gas Analyzers: Data Errors May Occur Medical Device Ongoing Action

Published: Monday, August 3, 2015
Last Updated: Monday, August 10, 2015

UMDNS Terms:

- Analyzers, Point-of-Care, Whole Blood, Gas/pH [18510]
- Analyzers, Point-of-Care, Whole Blood, Gas/pH/Electrolyte [18511]
- Analyzers, Point-of-Care, Whole Blood, Gas/pH/Electrolyte/Metabolite [18853]
- Analyzers, Point-of-Care, Whole Blood, Multianalyte [18509]

Product Identifier:

RAPIDPoint 500 Blood Gas Analyzers [Capital Equipment]
Siemens Material Nos.: 10492730, 10696855, 10696857, 10697306

Geographic Regions: Worldwide

Manufacturer(s): Siemens Healthcare Diagnostics Inc 511 Benedict Ave, Tarrytown, NY 10591, United States

Suggested Distribution: Clinical/Biomedical Engineering, Clinical Laboratory/Pathology, Pulmonology/Respiratory Therapy, Information Technology, Point-of-Care Coordination

Problem:

In a July 2015 Urgent Field Safety Notice letter submitted by an ECRI Institute member hospital, Siemens states when both ports (serial and Ethernet) are configured to transmit data, the data stream from one port may affect the data stream from the other port on the above analyzers. This problem may cause the message received by the laboratory information system (LIS) to include the following:

- Duplicate data,
- Missing data, or
- Data from a different patient (Siemens has not observed this scenario and the expected frequency is extremely unlikely).

Siemens states that this problem occurs only if the dual port LIS transmission feature is enabled on the above systems. If this problem occurs, only one analyte has the potential to be affected (i.e., if a panel of analytes is ordered, only one could be affected). Siemens states that a duplicate or missing analyte is of no risk to health. Siemens also states that a transposed sodium, potassium, glucose, lactate, carboxyhemoglobin, methemoglobin, or neonatal bilirubin result may lead to a potential risk to health when the true value is critical.

Action Needed:

Identify any affected product in your inventory. If you have affected product, verify that you have received the July 2015 Urgent Field Safety Notice letter and Field Correction Effectiveness Check Form from Siemens. Determine if your system is configured to use dual port LIS transmission. If your system is currently set up with dual LIS transmission enabled, one of the ports (either serial or Ethernet) must be disabled. To disable the connections, refer to the instructions and figures in the [letter](#). Complete the Field Correction Effectiveness Check form, and return it to Siemens using the instructions on the form. Siemens states that this problem will be resolved in the next software release. Keep a copy of the Urgent Field Safety Notice letter with your facility's laboratory records, and forward a copy of the letter to any facility to which you have further distributed affected product.

For Further Information:

Siemens customer care center

Website: [Click here](#)

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2015 Jul 30. Member Hospital. Siemens letter submitted by an ECRI Institute member hospital. (includes reply form) [Download](#)
- 2015 Aug 3. Manufacturer. The manufacturer confirmed the information in the source material
- 2015 Aug 10. MHRA FSN. 2015/007/030/601/001 [Download](#)
- 2015 Aug 10. MHRA FSN. 33114 Rev A. (includes reply form) [Download](#)

©2015 ECRI Institute
5200 Butler Pike, Plymouth Meeting, PA 19462-1298, USA
May be reproduced by subscribing institution for internal distribution only.