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ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Tuesday, January 13, 2015 Page 1

[High Priority] - A23471 01 : Zimmer—Trilogy Bone Screws: May Fracture during Insertion Medical Device Ongoing Action

UMDNS Terms:

- Screws, Bone [16101]

Product Identifier:

Trilogy Bone Screws:	Product Nos.:	Lot Nos.:
6.5 x 25 mm	00-6250-065-25	62784617, 62784618, 62784619, 62784621, 62784622, 62818701, 62818702, 62818703, 62818707, 62818709, 62818710, 62825611, 62836984
6.5 x 35 mm	00-6250-065-35	62748089, 62754329, 62793494, 62793495, 62793501, 62793502, 62793503, 62813612, 62813613

[Consumable]

Units distributed between August 19 and November 18, 2014

Geographic Regions: (Impact in additional regions has not been identified or ruled out at the time of this posting). Albania, Australia, Austria, Belgium, Brazil, Bulgaria, Canada, Canary Islands, Chile, China, Colombia, Czech Republic, Denmark, Egypt, Finland, France, Germany, Greece, Hungary, Ireland, Italy, Japan, Korea, Malaysia, The Netherlands, Norway, Poland, Portugal, Russia, Saudi Arabia, Serbia, Singapore, Slovakia, South Africa, Spain, Sweden, Switzerland, Taiwan, Turkey, U.K., U.S.

Manufacturer(s): Zimmer Inc 1800 W Center St, PO Box 708, Warsaw, IN 46581-0708, United States

Suggested Distribution: OR/Surgery, Orthopedics, Materials Management

Summary:

This alert provides additional information based on FDA Center for Devices and Radiological Health (CDRH) and U.K. Medicines and Healthcare Products Regulatory Agency (MHRA) source material regarding [Alert Accession No. A23471](#). FDA's CDRH states that the manufacturer initiated a recall by Urgent Medical Device Recall letter dated November 2014. Additional information is provided in the Geographic Regions field (see bolded region).

Problem:

[December 4, 2014]

In a November 21, 2014, Urgent Medical Device Recall Lot Specific letter submitted by an ECRI Institute member hospital, Zimmer states that the above screws may fracture during insertion. This problem may result in the following:

- If a screw fractures and does not protrude into the space occupied by an acetabular liner, a minor surgical delay may occur. If the fractured screw does protrude into the space occupied by an acetabular liner, a more lengthy delay may occur as the screw would have to be removed.
- If a screw fractures while being used to affix a knee femoral cut guide, a minor surgical delay may occur.
- If the fractured screw protrudes into the liner space, efforts to remove the screw and/or shell may result in pelvic bone damage, compromised fixation, and subsequent implant migration and revision.

The manufacturer has not confirmed the information provided in the source material.

Action Needed:

The following actions are those listed in [Alert Accession No. A23471](#). Identify, isolate, and discontinue use of any affected product in your inventory. If you have affected product, verify that you have received the November 21, 2014, Urgent Medical Device Recall Lot Specific letter from Zimmer. Your Zimmer representative will contact your facility to arrange for product removal. Notify all relevant personnel at your facility of the information in the letter, and forward a copy of the letter to any facility to which you have further distributed affected product. Zimmer recommends that you continue your normal post-operative follow-up routine for patients who have previously had affected product implanted. Inform Zimmer of any adverse events associated with use of affected product by e-mail at zimmer.per@zimmer.com. U.S. customers should also report any adverse events to FDA's MedWatch Adverse Event Reporting program by fax at (800) 332-0178 ; by mail (using postage-paid FDA Form 3500, available [here](#)) at Food and Drug Administration, 5600 Fishers Lane, Rockville, MD 20852-9787; or online at the [MedWatch website](#).

For Further Information:

Zimmer customer call center

Tel.: (800) 348-2759, 8 a.m. to 8 p.m. Eastern time Monday through Friday, or (877) 946-2761

Website: [Click here](#)

References:

- United States. Food and Drug Administration. Center for Devices and Radiological Health. Class 2 device recall Trilogy bone screws [online]. 2014 Dec 17 [cited 2015 Jan 2]. Available from Internet: [Click here](#).

©2015 ECRI Institute
5200 Butler Pike, Plymouth Meeting, PA 19462-1298, USA
May be reproduced by subscribing institution for internal distribution only.

- Great Britain. Medicines and Healthcare Products Regulatory Agency. Zimmer. Osteosynthesis. Trilogy bone screws—6.5mm X 35mm & 6.5mm X 25mm [online]. London: Department of Health; 2014 Dec 22 [cited 2015 Jan 5]. (Field safety notice; reference no. 2014/012/015/081/002). Available from Internet: [Click here](#) .

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#) .

Source(s):

- 2014 Dec 4. Member Hospital. Zimmer letter submitted by an ECRI Institute member hospital. [Download](#)
- 2015 Jan 2. FDA CDRH Database. Class II. Z-0824-2015 [Download](#)
- 2015 Jan 5. MHRA FSN. 2014/012/015/081/002 [Download](#)
- 2015 Jan 5. MHRA FSN. (includes reply form) [Download](#)