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## ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Friday, June 12, 2015 Page 1

### [High Priority] - A24431 : Intuitive—5 mm Thoracic Graspers Used with da Vinci S and da Vinci Si Surgical Systems: Contact/Abrasive Forces May Cause Damage to Main Tube Coating Medical Device Ongoing Action

Published: Thursday, May 28, 2015

#### UMDNS Terms:

- Telesurgery Systems, Surgical, Minimally Invasive [18600]

#### Product Identifier:

5 mm Thoracic Graspers used with the following da Vinci Surgical Systems: (1) S, (2), Si [Capital Equipment]  
Thoracic Grasper Part Nos.: 420343-01, 420343-02

**Geographic Regions:** (Impact in additional regions has not been identified or ruled out at the time of this posting), Belgium, Brazil, Canada, Chile, Colombia, Denmark, France, Germany, India, The Netherlands, Panama, Singapore, South Korea, Switzerland, Taiwan, Turkey, U.K., U.S.

**Manufacturer(s):** Intuitive Surgical Inc 1266 Kifer Rd Bldg 101, Sunnyvale, CA 94086-5304, United States

**Suggested Distribution:** Cardiology/Cardiac Catheterization Laboratory, Clinical/Biomedical Engineering, OR/Surgery, Risk Management/Continuous Quality Improvement, Materials Management

#### Problem:

In a May 21, 2015, Field Safety Notice Urgent Medical Device Recall letter submitted by ECRI Institute member hospitals, Intuitive states that contact/abrasive forces applied to the main tube coating of the above thoracic graspers may cause the damage to the main tube (for an image of damage to the main tube see Figure 1 in the [Field Safety notice Urgent Medical Device Recall letter](#)). Intuitive also states that the force being applied can come from a damaged (e.g., bent, abraded, notched) cannula or from contact with an adjacent instrument. Intuitive further states that if the main tube is damaged during a procedure, fragments of the main tube coating may fall into the patient; however, the user may visually detect instrument failure either through the surgeon console or at the patient side during instrument exchanges. If damage is observed, the user would need to replace the instruments, resulting in a minor delay in the procedure, and material that fell into the patient would need to be retrieved or removed through suction-irrigation. Intuitive states that the material that may fragment from the main tube of the above thoracic graspers meets the ISO-10993 biocompatibility and implantability testing requirements and manifests no signs of toxicity. Intuitive also states that for thoracic and extracardiac surgeries, main tube damage with fragments of coating material falling into the surgeon's view could have a minor immediate health effect and a temporary delay in the procedure could occur while the fragment is retrieved and the instrument replaced. For intracardiac surgery, if a fragment falls into the patient during the procedure, it is likely to be seen because of the smaller operating space and the high contrast between the dark coating material and the light intracardiac tissue and would be easily retrieved and the instrument replaced with only a brief delay in the procedure. If a particle becomes enclosed within the heart, the patient could be at risk for an embolic event such as a stroke. Intuitive further states that this high risk is mitigated by the Cardiac Surgical Standard Operating Procedure of inspecting the atrium and valve before closing to ensure that no material of any kind is left behind. Intuitive states that it has received no reports of patient injuries as a result of this problem. Intuitive also states that this problem affects only the above thoracic graspers; thoracic graspers with part number 420343-03 are not affected. The manufacturer has not confirmed the information provided in the source material.

#### Action Needed:

Identify any affected product in your inventory. For the location of the part and version number on the instrument, see Figure 2 in the [Field Safety Notice Urgent Medical Device Recall letter](#). If you have affected product, verify that you have received the May 21, 2015, Field Safety Notice Urgent Medical Device Recall letter and Acknowledgment Form from Intuitive. Notify all relevant personnel at your facility of the information in the letter. To arrange for product return and to obtain return material authorizations (RMAs) contact the Intuitive customer service department using the information below. Upon receipt of affected product, Intuitive will provide your facility with new replacement thoracic graspers (part number 420343-03). To determine the number of replacement instruments to ship to your facility, Intuitive will determine the total combined number of remaining lives of all instruments returned and round the number of lives up to the nearest increment of 20 (i.e., if a customer returned 3 instruments with 15 lives each, the total number of remaining lives is 45. The 45 remaining lives would be rounded up to 60 lives, and 3 thoracic instruments, each with 20 lives, would be shipped as replacement instruments). Complete the Acknowledgment Form, and return it to Intuitive using the instructions on the form. Retain a copy of the letter with your records.

#### For Further Information:

Intuitive local representative or Intuitive Surgical customer service department  
**North America and South America:**  
Tel.: (800) 876-1310 (select option 3), 6 a.m. to 5 p.m. Pacific time  
**Japan:**  
Tel.: (0120) 565635 or (03) 55751362, 9 a.m. to 6 p.m. Japan time  
**South Korea:**  
Tel.: (02) 32713200, 9 a.m. to 6 p.m. Korea time  
**Europe, Middle East, Africa, and Asia:**  
Tel.: (800) 08212020 or 41 (21) 8212020, 8 a.m. to 6 p.m. Central European time  
E-mail: [ics@intusurg.com](mailto:ics@intusurg.com)  
Website: [Click here](#)

#### Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

#### Source(s):

- 2015 May 27. Member Hospital. Intuitive letter submitted by ECRI Institute member hospitals: 2955842-05-07-2015-005-R (includes reply form) [Download](#)

©2015 ECRI Institute  
5200 Butler Pike, Plymouth Meeting, PA 19462-1298, USA  
May be reproduced by subscribing institution for internal distribution only.