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## ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Monday, June 1, 2015 Page 1

### [Normal Priority ] - S0274 : Stryker—InTouch Beds: Pendant Cord May Be Damaged if Routed Between Side Rail and Bed Frame [ECRI Institute User Experience Network] Medical Device Special Report

Published: Thursday, May 14, 2015

**UMDNS Terms:**

- Beds, Electric [10347]

**Product Identifier:**

InTouch Critical Care Beds [Capital Equipment]

**Geographic Regions:** Worldwide

**Manufacturer(s):** Stryker Corp 2825 Airview Blvd, Kalamazoo, MI 49002, United States

**Suggested Distribution:** Clinical/Biomedical Engineering, Critical Care, Nursing

**Problem:**

- If the pendant cord (which connects the hand pendant control to the bed frame) is routed between the side rail and the frame of the above beds, the cord insulation may fray due to the side rail rubbing against the cord during repeated raising and lowering of the rails.
- The cord may become damaged to the point where wires will be exposed.

**Manufacturer's Corrective Action/Recommendations:**

- Stryker states that the solution to this problem is typically to route the cord differently. Stryker does not install the pendant/cord, but can assist by sending a service technician to the hospital to help find better cord routes.
- The pendant is a low voltage pendant so there is no danger to touching exposed wires on the pendant cord. Stryker has also stated that the pendant ceasing to function properly will not affect other bed functions.

**ECRI Institute Recommendations:**

1. Do not route the pendant cord between the bed side rails and frame, and alert all users of the Stryker InTouch Critical Care Bed that the pendant cord insulation may fray and expose the cord wires if it is routed between the bed side rails and frame.
2. If you are having issues finding a proper route for the pendant cord, contact Stryker Medical Technical Support by telephone at (800) 327-0770.

**Background:**

- ECRI Institute received multiple reports of Stryker InTouch Critical Care Bed pendants being damaged by being caught in between the side rail and frame.
- One of the reporting hospitals now routes the cord under and through a gap in the side rail (Figure 1), which has alleviated the problem. However, with this new route, the cord occasionally snags door handles or other objects when transporting the bed.

*Figure 1. Updated cord routing strategy used by reporting hospital.*

**Comments:**

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

**Source(s):**

- 2015 May 14. ECRI Institute researched member report.
- 2015 May 14. Figure 1. [Download](#)

©2015 ECRI Institute  
5200 Butler Pike, Plymouth Meeting, PA 19462-1298, USA  
May be reproduced by subscribing institution for internal distribution only.