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ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Wednesday, December 14, 2016 Page 1

[High Priority] - A27550 : NIHON KOHDEN—cardiolife AED-2100K Automated External Defibrillators: Battery May Become Depleted Prematurely after Post-Self-Test Shutdown Failure Medical Device Ongoing Action

Published: Wednesday, December 7, 2016

UMDNS Terms:

- Defibrillators, External, Automated [17116]

Product Identifier:

cardiolife AED-2100K Automated External Defibrillators (AEDs) [Capital Equipment]

Geographic Regions: Asia, Central America, Europe, Middle East, South America

Manufacturer(s): Nihon Kohden Europe GmbH Raiffeisenstrasse 10, D-61191 Rosbach vor der Hoehe, Germany

Suggested Distribution: Cardiology/Cardiac Catheterization Laboratory, Clinical/Biomedical Engineering, Critical Care, Emergency/Outpatient Services, Nursing, Information Technology, Home Care, EMS/Transport

Problem:

In a December 2016 Important Field Safety Notice letter posted by the German Federal Institute for Drugs and Medical Devices (BfArM), NIHON KOHDEN states that the power on the above AEDs may not turn off after completion of the daily internal self-test, potentially resulting in the battery becoming discharged within one day. NIHON KOHDEN also states that if this problem occurs, the status indicator remains lit (red) but no alarm sound is generated and defibrillation may not be available when required because of a depleted battery.

Action Needed:

Identify any affected AEDs in your inventory. If you have affected AEDs, verify that you have received the December 2016 Important Field Safety Notice letter, receipt of Field Safety Notice form, and Return Delivery Sheet from NIHON KOHDEN. Notify all relevant personnel at your facility of the information in the Important Field Safety Notice letter. Complete the receipt of the Field Safety Notice form, and return it to NIHON KOHDEN using the information in the letter. Upon receipt of the form, a NIHON KOHDEN representative will contact your facility to arrange to pick up affected AEDs. Complete the Return Delivery Sheet, and return it, along with affected product, to NIHON KOHDEN. NIHON KOHDEN will update the software to version 02-02.

For Further Information:

NIHON KOHDEN technical department

Tel.: 49 (6003) 827160

E-mail: NKE-SERVICE2@nke.de

Website: [Click here](#)

References:

- Germany, Federal Institute for Drugs and Medical Devices. Urgent field safety notice for automated external defibrillators cardiolife AED-2100K, NIHON KOHDEN EUROPE GmbH [online]. 2016 Nov 8 [cited 2016 Nov 28]. Available from Internet: [Click here](#).

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2016 Nov 28. BfArM (Germany). 09428/16 [Download](#)
- 2016 Nov 28. BfArM (Germany). Nihon Kohden Reference No. AED-2100K/2016/1 [Download](#)
- 2016 Dec 7. Manufacturer. The manufacturer confirmed the information in the source material.