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ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Monday, November 21, 2016 Page 1

[High Priority] - A27566 : GE— Avance and Amingo Anesthesia Devices with Large Tray Insert Accessory: May Transition to System Malfunction State If Storage Drawer Is Closed Forcefully Medical Device Ongoing Action

Published: Monday, November 14, 2016

UMDNS Terms:

- Anesthesia Units [10134]

Product Identifier:

Anesthesia Devices with Large Tray Insert Accessory: (1) Amingo, (2) Avance, (3) Avance CS² [Capital Equipment]
Large Tray Insert Part No. 1009-3260-000

Geographic Regions: Worldwide

Manufacturer(s): GE Healthcare9900 Innovation Dr, Wauwatosa, WI 53226, United States

Suggested Distribution: Anesthesia, Clinical/Biomedical Engineering, OR/Surgery, Pulmonology/Respiratory Therapy

Problem:

In a November 3, 2016, Urgent Medical Device Correction letter submitted by ECRI Institute member hospitals, GE states that the above devices may transition to a system malfunction state if the lower storage drawer containing the optional large tray insert accessory is closed with an abnormally high amount of force, potentially resulting in a loss of patient ventilation and hypoxia if the system malfunction is left unresolved. GE also states that the above devices will operate as follows in the system malfunction state:

- Alternative oxygen flow will automatically activate within a few seconds.
- High-priority audio and visual alarms will activate.
- Instructions to set the oxygen flow and manually ventilate the patient will display.
- Anesthetic agent will continue to be delivered at the existing vaporizer setting.

GE further states that it has received no reports of injury as a result of this problem.

Action Needed:

Identify any affected systems in your inventory. If you have affected systems, verify that you have received the November 3, 2016, Urgent Medical Device Correction letter from GE. Remove the optional large tray insert accessory if it is installed. You can continue using affected devices once the large tray insert accessory is removed. GE recommends destroying all tray insert accessories. You may choose to return the tray insert to GE. To return the tray insert to GE, e-mail a request to Avance_Drawer.Insert_RMA@ge.com.

For Further Information:

GE
Tel.: (800) 437-1171
Website: [Click here](#)

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2016 Nov 10. Member Hospital. GE Reference No. 34079 [Download](#)
- 2016 Nov 14. Manufacturer. The manufacturer confirmed the information provided in the source material.

©2016 ECRI Institute
5200 Butler Pike, Plymouth Meeting, PA 19462-1298, USA
May be reproduced by subscribing institution for internal distribution only.