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## ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Monday, November 21, 2016 Page 1

### [Normal Priority ] - H0348 : Intuitive—da Vinci EndoWrist Vessel Sealers: Keys to Preventing Malfunctions and Patient Harm [ECRI Exclusive Hazard Report] Medical Device Hazard Report

Published: Thursday, November 17, 2016

#### UMDNS Terms:

- Telesurgery Systems, Surgical, Minimally Invasive [18600]

#### Product Identifier:

da Vinci EndoWrist Vessel Sealers used with the following da Vinci Surgical Systems: (1) Si, (2) Si-e, (3) Xi [Consumable, Capital Equipment]  
Vessel Sealer Part Nos.: 410322, 480322

Geographic Regions: Worldwide

Manufacturer(s): Intuitive Surgical Inc 1020 Kifer Rd, Sunnyvale, CA 94086-5304, United States

Suggested Distribution: Clinical/Biomedical Engineering, OR/Surgery, Risk Management/Continuous Quality Improvement

#### Problem:

1. Incorrect use of the da Vinci EndoWrist Vessel Sealer instrument can cause the following:
  - Incomplete cutting or sealing of tissue
  - Damage to vessel sealer jaws
  - Detachment of the knife cable from knife blade
2. These problems can cause the following device malfunctions:
  - The vessel sealer may not fully seal tissue, despite audible tones indicating the end of the sealing cycle.
  - The vessel sealer blades may not fully transect tissue or the blade may not fully retract and the system may display a blade exposure error.
  - The blades may detach from the instrument.
3. Some complications may not be immediately detectable and can cause patient harm, such as the following:
  - Insufficient vessel sealing may lead to unexpected bleeding that may not be readily apparent.
  - The exposed blades may cause inadvertent patient injury if the user continues to use the device with the blade exposed and the blade contacts tissue.

#### Manufacturer's Recommendations:

1. Only use the da Vinci EndoWrist Vessel Sealer instrument to cut vessels up to 7 mm and tissue bundles that fit within the jaws. Skeletonize large tissue bundles whenever possible.
2. Before activating the knife blade, ensure that the vessel sealer jaws are no more than 4 mm apart.
3. Do not use the vessel sealer on hard objects such as bone, calcified tissue, staples, or clips.
4. Do not cut sealed tissue unless the desired tissue effect is observed. Blanching, smoke, and/or bubbling are all indications of energy delivery.
5. Follow on-screen prompts for resolving blade exposure: if prompted, remove and re-install the vessel sealer. If the blade exposure resolves upon insertion, use may be continued.
6. See Intuitive's [Updated Customer Letter](#) .

#### ECRI Institute Recommendations:

While ECRI Institute concurs with all of the above recommendations, we have not examined any affected instruments. In addition to following the manufacturer's recommendations listed above, ECRI Institute recommends the following:

1. Monitor any changes in sealing time as this may indicate device failure.

©2016 ECRI Institute  
5200 Butler Pike, Plymouth Meeting, PA 19462-1298, USA  
May be reproduced by subscribing institution for internal distribution only.

2. Inspect tissue for adequate sealing.
3. Do not continue to use an instrument if malfunction is suspected.

**Background:**

- da Vinci EndoWrist Vessel Sealer instruments are used with da Vinci Si, da Vinci Si-e, and da Vinci Xi systems to cut and cauterize tissue.
- The instruments include independent cut and seal functions.
- Cutting is achieved using a cable-driven knife blade that extends upon activation and retracts once cutting is complete.

**Comments:**

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

**Source(s):**

- 2016 Nov 10. ECRI Institute
- 2016 Nov 11. Manufacturer Letter. Document submitted by the manufacturer [Download](#)