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ANEXO

www.ecri.org . Printed from *Health Devices Alerts* on Tuesday, August 1, 2017 Page 1

[High Priority] - A28950 : Otto Bock—Voyager Evo Wheelchairs: Caster Fork May Turn Out of Position

Medical Device Ongoing Action

Published: Monday, July 24, 2017

Last Updated: Thursday, July 27, 2017

UMDNS Terms:

- Wheelchairs, Mechanical [16949]

Product Identifier:

Voyager Evo Wheelchairs [Capital Equipment]

Product No. 480A76=00000_K

Geographic Regions: Argentina, Australia, Austria, Belgium, Canada, China, Colombia, Czech Republic, Finland, France, Germany, Hong Kong, Hungary, Ireland, Italy, Korea, Luxembourg, Mexico, The Netherlands, Poland, Romania, Russia, Singapore, Slovakia, South Africa, Switzerland, Thailand, U.K.

Manufacturer(s): Otto Bock HealthCare GmbH Max-Naeder-Strasse 15, Duderstadt, D-37115, Germany

Suggested Distribution: Clinical/Biomedical Engineering, Facilities/Building Management, Home Care, Physical Therapy/Rehabilitation

Problem:

In a July 7, 2017, Urgent Safety Information letter posted by the German Federal Institute for Drugs and Medical Devices (BfArM), Otto Bock states that the caster forks may turn out of position on the above wheelchairs, potentially causing connections to become loose and the castor fork to go out of adjustment. This may negatively influence driving characteristics and potentially lead to the front of the wheelchair dropping down and falling forward. Otto Bock also states that this problem may be caused by adjusting the caster journal angle after optimizing the user's sitting position.

Action Needed: Identify any affected wheelchairs in your inventory and identify any patients using affected wheelchairs. If you have affected wheelchairs and/or affected patients, verify that you have received the July 7, 2017, Urgent Safety Information letter from Otto Bock. Regularly check the alignment of the caster fork axis and watch for changed driving characteristics and/or play between the components during use. If you identify this failure mode on a wheelchair, contact your Otto Bock local distributor. Inform all users of the information in the letter.

For Further Information:

Otto Bock customer care center

Tel.: 49 (69) 9999394

E-mail: ccc@ottobock.de

Website: [Click here](#)

References:

- Germany. Federal Institute for Drugs and Medical Devices. Urgent field safety notice for Voyager Evo Starrrahmenrollstuhl by Otto Bock Mobility Solutions [online]. 2017 Jul 14 [cited 2017 Jul 19]. Available from Internet: [Click here](#).

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2017 Jul 19. BfArM (Germany). 06599/17 [Download](#)
- 2017 Jul 19. BfArM (Germany). [Download](#)
- 2017 Jul 24. Manufacturer. Manufacturer confirmed information