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ANEXO

www.ecri.org . Printed from *Health Devices Alerts* on Friday, August 18, 2017 Page 1

[High Priority] - H0394 : ICU Medical/Q Core Medical—Sapphire Lockboxes: Can Be Defeated, Providing Medication Access to Unauthorized Persons [ECRI Exclusive Hazard Report]

Medical Device Hazard Report

Published: Tuesday, August 15, 2017

UMDNS Terms:

- Infusion Pumps, Analgesic, Patient-Controlled [16924]

Product Identifier:

Sapphire Lockboxes [Capital Equipment]

Geographic Regions: Worldwide

Distributor(s): • ICU Medical Inc951 Calle Amanecer, San Clemente, CA 92673, United States

Manufacturer(s): Q Core Medical Ltd29 Yad Haruzim Street, Netanya IL-42505, Israel

Suggested Distribution: Anesthesia, Clinical/Biomedical Engineering, Critical Care, Nursing, Risk Management/Continuous Quality Improvement, Pharmacy, IV Therapy

Problem:

1. ICU Medical/Q Core Medical Sapphire lockboxes may be defeated with items available in a care unit.
2. Patients or other unauthorized persons (e.g., family members) may tamper with the lockbox without breaking it and gain access to medication.
 - The ability to gain access to the medication presents a risk of unauthorized persons tampering with and/or stealing the medication.
3. Interruption and/or delay in pain therapy may occur as a result of this problem.
4. The pump does not recognize or alarm if the lockbox is tampered with, even if an infusion is ongoing.

ECRI Institute Recommendations:

1. Notify clinician staff about the potential for unauthorized persons to gain access to medication when an ICU Medical/Q Core Medical lockbox is used.
2. Clinicians should conduct a physical and visual check at a specific interval (e.g., every shift change) to ensure that the lockbox is appropriately secured.
3. If unauthorized access to medication is suspected, report the incident per your facility's policy.

Manufacturer's Perspectives or Comments:

1. The ICU Medical/Q Core Medical operational manual states:
 - "PCA Lockboxes are designed to secure the IV bag, primarily for treatments involving narcotics or opioids, without interrupting the treatment workflow."
2. There are no known cases of medication diversion when locked in the Sapphire lockbox.

Background:

- ECRI Institute has received reports from member facilities about patients gaining access to patient-controlled analgesia (PCA) pump lockboxes.
 - Reports pertain to a different manufacturer's PCA pump lockbox than what is discussed in this report.
- ECRI Institute identified the ability to gain access to the ICU Medical/Q Core Medical lockbox without a key in our laboratory testing.
- ECRI Institute is currently not aware of any alternative methods to secure the medication, beyond the lockbox, with this particular pump model.
- In the past, ECRI Institute has mainly received reports and investigated incidents regarding medication diversion by hospital staff. Recent reports of patients gaining access to their medications suggest a growing trend. Considering the ongoing opioid epidemic in the U.S., we want to ensure that facilities are aware of these shortcomings and provide them with recommendations to act accordingly.

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or

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www.ecri.org . Printed from *Health Devices Alerts* on Friday, August 18, 2017 Page 2

source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#) .

Source(s):

- 2017 Aug 1. ECRI Institute researched member report