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ANEXO

www.ecri.org . Printed from *Health Devices Alerts* on Tuesday, September 5, 2017 Page 1

[High Priority] - A29086 : Fujifilm—Various FDR/CR Digital Mammography Systems: May Assign Duplicated ID Number to Exposure Images, Potentially Causing Previous Studies to Be Overwritten Medical Device Ongoing Action

Published: Wednesday, August 16, 2017
Last Updated: Thursday, August 17, 2017

UMDNS Terms:

- Radiographic Systems, Digital, Mammographic [18432]

Product Identifier:

Digital Mammography Systems: (1) FDR-1000AWS; (2) FDR-2000AWS; (3) FDR-3000AWS; (4) CR-IR363AWS [Capital Equipment]
Software Versions: V5.0, V5.1, V5.2, V6.0, V6.1, V7.0
Product Nos.: 00000000006730135, 16830333, 16830335, 16940041, 16940043, 46620199, 56320374, 56320375, 56320377, 56520409, 56720475, 56720476, 57120599, 57220619, 66420755, 66520794, 66520795, 66620829, 66620830, 66620831

Geographic Regions: (Impact in specific regions has not been identified or ruled out at the time of this posting), Worldwide

Manufacturer(s): FUJIFILM UK Ltd Unit 10A St Martins Business Centre, Bedford, MK42 0LF, England

Suggested Distribution: Clinical/Biomedical Engineering, Obstetrics/Gynecology/Labor and Delivery, Oncology, Diagnostic Imaging

Problem:

In a Field Safety Notice letter posted by the U.K. Medicines and Healthcare Products Regulatory Agency (MHRA), Fujifilm states that the above software versions for the above digital mammography systems may assign a duplicated ID number to exposure images. If an image with this error is transmitted to PACS, it may overwrite the image already stored on PACS. Fujifilm also states that this problem may cause an unexpected image to appear when referring to a previous study; however, the image and data of a new study are unaffected, including a diagnosis. Fujifilm further states that depending on PACS specifications, images may not be overwritten or a transmission error will occur. The manufacturer has not confirmed the information provided in the source material.

Action Needed:

Identify any affected systems in your inventory. If you have affected systems, verify that you have received the Field Safety Notice letter and Customer Feedback Form from Fujifilm. Ensure that the latest images transmitted to PACS exist in the new study; if the image with a duplicated ID number assigned is transmitted to PACS, it will not be found in a newly created study. The image transmitted to PACS with this error overwrites an existing (older) study and is stored in another (older) study. If this error occurs, contact your local Fujifilm office. Fujifilm will contact your facility to arrange for correction of this problem.

For Further Information:

Fujifilm
Website: [Click here](#)
References:

Great Britain. Medicines and Healthcare Products Regulatory Agency. FUJIFILM: FDR-1000AWS/FDR-2000AWS/FDR-3000AWS/CR-IR363 [online]. London: Department of Health; 2017 Aug 14 [cited 2017 Aug 16]. (Field safety notice; reference no. 2017/008/010/701/020. Available from Internet: [Click here](#) .

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#) .

Source(s):

- 2017 Aug 16. MHRA FSN. 2017/008/010/071/020 [Download](#)
- 2017 Aug 16. MHRA FSN. FFID_20170707 (includes reply form) [Download](#)