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### ANEXO

www.ecri.org . Printed from *Health Devices Alerts* on Wednesday, August 30, 2017 Page 1

**[High Priority ] - A29133 : Olympus—Cystoscopy Bridges and Working Inserts: Adhesive Fragments May Detach**  
**Medical Device Ongoing Action**

Published: Thursday, August 24, 2017

**UMDNS Terms:**

- Cystoscopes [11112]
- Surgical Retractor Kits, Urinary Tract [33236]

**Product Identifier:**

Products:	Model Nos.:	Lot Nos.:
Bridge, One Way	A20976A	146W through 172W
Bridge, Two Way	A20977A	146W through 174W
Working Insert with Ramp, One Way	A20975A	13ZW through 172W

[Consumable]

**Geographic Regions:** Worldwide

**Manufacturer(s):** Olympus Winter & Ibe GmbHKuehnstrasse 61, Hamburg, D-22045, Germany

**Suggested Distribution:** OR/Surgery, Urology, Materials Management

**Problem:** In an August 18, 2017, Urgent Medical Device Removal and Corrective Action letter submitted by an ECRI Institute member hospital, Olympus states that it has received reports of incidents in which fragments of adhesive detached from inside the working channel of the above bridges. Olympus also states that it has received reports of cracking, chipping, missing pieces, and delamination of the adhesive. If adhesive detaches during the intended use of the above bridges or working inserts (e.g., when inserting an instrument through the working channel), a fragment of adhesive may fall inside the patient's bladder or urethra, necessitating retrieval. Although they are typically flushed out with irrigation fluid or passed naturally, the retrieval of large fragments of adhesive could prolong the procedure or require additional surgical treatment. Olympus further states that it has received no reports of patient injury related to this problem.

**Action Needed:**

Identify and discontinue use of any affected product in your inventory. If you have affected product, verify that you have received the August 18, 2017, Urgent Medical Device Removal and Corrective Action letter and Reply Form from Olympus. To arrange for product return, bridge rework, and working insert replacement, contact the Olympus customer solutions department using the information below. Olympus states that, based on its investigation results and risk assessment, affected product is safe to use until reworked or replaced. Complete the Reply Form and return it to Olympus using the instructions on the form. Forward a copy of the letter to any facility to which you have further distributed affected product, and notify Olympus of the transfer.

**For Further Information:**

Olympus customer solutions department  
 Tel.: (800) 848- 9024 (select option 3)  
 Website: [Click here](#)

**Comments:**

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

**Source(s):**

- 2017 Aug 23. Member Hospital. Olympus letter submitted by ECRI Institute member hospital (includes reply form) [Download](#)
- 2017 Aug 24. Manufacturer. Manufacturer confirmed information