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### ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Friday, June 23, 2017 Page 1

#### [High Priority ] - A28765 : ArjoHuntleigh—Model ACS900 Flowtron Pumps: May Deliver Constant Pressure to Calf Garment without Triggering Alarm Medical Device Ongoing Action

Published: Thursday, June 22, 2017

#### UMDNS Terms:

- Circulatory Assist Units, Peripheral Compression, Intermittent [10969]

**Product Identifier:** Model ACS900 Flowtron Pumps [*Capital Equipment*]  
Software Version V1.099; Serial Nos.: 1400028031 through 1600048470 Units manufactured between September 26, 2014, and December 20, 2016

**Geographic Regions:** Australia, Austria, Bahrain, Brazil, Canada, Chile, Colombia, Cyprus, Czech Republic, Denmark, Ecuador, France, Germany, Greece, Hong Kong, Hungary, India, Ireland, Israel, Italy, Japan, Jordan, Kuwait, Lebanon, Malaysia, Mauritius, Mexico, The Netherlands, New Zealand, Norway, Oman, Pakistan, Paraguay, Poland, Qatar, Saudi Arabia, Singapore, South Africa, Spain, Sweden, Switzerland, Tanzania, United Arab Emirates, U.K., U.S.

**Manufacturer(s):** ArjoHuntleigh Polska Sp z o o ulica Ks Wawrzyniaka 2, PL-62-052 Komorniki, Poland

**Suggested Distribution:** Cardiology/Cardiac Catheterization Laboratory, Clinical/Biomedical Engineering, Critical Care, Emergency/Outpatient Services, Nursing, OR/Surgery, Information Technology, Home Care

#### Problem:

In an Urgent Field Safety Notice letter posted by the U.K. Medicines and Healthcare Products Regulatory Agency (MHRA), ArjoHuntleigh states that the above pumps may deliver constant pressure to the garment without triggering the alarm, potentially resulting in patient injury. ArjoHuntleigh also states that it has received reports of this problem occurring in conjunction with a blank LCD display.

#### Action Needed:

Identify any affected pumps in your inventory. If you have affected pumps, verify that you have received the Urgent Field Safety Notice letter and Customer Response Form from ArjoHuntleigh. Upgrading affected pumps to software revision V2.000 will ensure more resilient systems. The firm states that you may continue to use affected pumps until an upgrade can be performed, on the condition that the patient's limbs are monitored frequently and the garments are checked to ensure that they are correctly fitted to the patient and that they are deflating on a regular bases. ArjoHuntleigh also recommends that you check the pump LCD display regularly to ensure that it indicates that the garments are inflating and deflating correctly. Notify all relevant personnel at your facility of the information in the Urgent Field Safety Notice letter, and forward a copy of the letter to any facility to which you may have further transferred affected product. Complete the Customer Response Form, and return it to ArjoHuntleigh using the information on the form. Upon receipt of the form, an ArjoHuntleigh representative will contact your facility to arrange for service and upgrade at no cost. Retain a copy of the letter with your records.

#### For Further Information:

ArjoHuntleigh  
Website: [Click here](#)

#### References:

- Great Britain. Medicines and Healthcare Products Regulatory Agency. ArjoHuntleigh: Flowtron ACS900 [online]. London: Department of Health; 2017 Jun 12 [cited 2017 Jun 19] (Field safety notice; reference no. 2017/006/006/701/008). Available from Internet: [Click here](#) .

#### Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#) .

#### Source(s):

- 2017 Jun 19. MHRA FSN. 2017/006/006/701/008 [Download](#)
- 2017 Jun 19. MHRA FSN. ArjoHuntleigh Reference No. FSN/SUZ/001-2017 (includes reply form) [Download](#)
- 2017 Jun 21. Manufacturer. The manufacturer confirmed the information in the source material.
- 2017 Jun 22. Health Canada Recall Listings. Type II. RA-63684 [Download](#)

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