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ANEXO 1

www.ecri.org . Printed from Health Devices Alerts on Monday, February 27, 2017 Page 1

[High Priority] - A27970 01 : Roche—cobas p 612 Preanalytical Systems: Manufacturer Revises Software to Correct Aspiration Process Medical Device Ongoing Action

Published: Wednesday, February 22, 2017

UMDNS Terms:

- Analyzers, Laboratory, Clinical Chemistry/Immunoassay [20821]

Product Identifier: cobas p 612 Preanalytical Systems [Capital Equipment]
GMDI/Part No. 07563116001
Serial Nos.: 63001216 and higher
Software Versions: 2.0.0, 2.0.1, 2.0.2, 2.0.3

Geographic Regions: Asia, Australia, Brazil, Chile, Colombia, Cuba, Ecuador, Europe, Jamaica, Panama, Philippines, Tunisia, U.K., U.S.

Manufacturer(s): Roche Diagnostics GmbH Sandhoferstrasse 116, D-68305 Mannheim, Germany

Suggested Distribution: Clinical/Biomedical Engineering, Clinical Laboratory/Pathology, Information Technology

Summary:

This Alert provides new information based on manufacturer correspondence and a December 15, 2016, Urgent Field Safety Notice letter posted by the U.K. Medicines and Healthcare Products Regulatory Agency (MHRA) regarding [Alert Accession No. A27970](#). New information is provided in the following fields:

- Product Identifier (see bolded information)
- Geographic Distribution (see bolded regions)
- Problem
- Action Needed

Problem:

[February 22, 2017]

In a December 15, 2016, Urgent Field Safety Notice letter posted by MHRA, Roche states that when the above system is using disposable non-filter tips, the system software may incorrectly calculate the submersion depth. An incorrect submersion depth calculated by the software when non-filter tips are used may lead to the aspiration of an increased number of aspirated clots, caused by the tip moving past the serum layer and entering the gel or blood cake layer. The clots are still being correctly detected, creating a system alarm. When a clot is detected, the pipettor may fail to stop the aspiration. This creates an underpressure in the pipettor tip, which may lead to liquid splashing up inside the pipettor tip when the blockage is removed, potentially resulting in a droplet forming at the pipettor nozzle. The firm received 1 report regarding an increased number of aspirated clots. While no reports of sample carryover have been received, the effect could not be excluded. The effect of cross-contamination on patient results depends on the assay (e.g., on the sensitivity of the assay) as well as the size of the droplet, the aspirated amount of sample, and the concentration of the primary sample from which the droplet originates. A medical risk cannot be excluded. The firm adds that all cobas p 612 (63x) systems that use disposable filter tips are not affected by a potential nozzle contamination or sample carryover. The submersion depth with filter tips is correct and, in addition, filter tips physically prevent liquid reaching and contaminating the nozzle.

[February 1, 2017]

FDA's Center for Devices and Radiological Health (CDRH) states that Roche is revising the software of the above systems to correct the aspiration process. FDA's CDRH also states that Roche initiated a product correction of the above systems by Urgent Medical Device Correction letter sent by fax or e-mail on December 27, 2016.

Action Needed: Identify any affected product in your inventory. If you have affected product, verify that you have received the December 15, 2016, Urgent Field Safety Notice letter from Roche. The manufacturer states that a new software version is available and that the upgrade of affected systems is planned for July 31, 2017. Until the new software version 2.0.4 is installed, stop aliquoting sample tubes if you are using affected systems for aliquoting in combination with non-filter tips. As a workaround, if sample aliquoting is required until the software version 2.0.4 is installed, systems using non-filter tips can be modified free of charge by your Roche local field service representative to use disposable filter tips. Inform all relevant personnel at your facility of the information in the letter, and forward a copy of the letter to any facility to which you have further distributed affected product.

For Further Information:

Roche technical support hotline
 Tel: (0808) 1001920 (U.K.), (1800) 409564 (Ireland)
 Website: [Click here](#)

References:

- United States. Food and Drug Administration. Center for Devices and Radiological Health. Class 2 device recall cobas p 612 preanalytical system [online]. 2017 Jan 24 [cited 2017 Feb 1]. Available from Internet: [Click here](#).
- Great Britain. Medicines and Healthcare Products Regulatory Agency. Roche: cobas p 612 pre-analytical system [online]. London: Department of Health; 2017 Feb 20 [cited 2017 Feb 22]. (Field safety notice; reference no. 2016/012/020/701/016). Available from Internet: [Click here](#).

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [FDA Format Guide](#).

Source(s):

- 2017 Feb 22. FDA CDRH Database. Class II. Z-1090-2017 [Download](#)
- 2017 Feb 22. Manufacturer. Manufacturer confirmed information contained in source material

©2017 ECRI Institute
 5200 Butler Pike, Plymouth Meeting, PA 19462-1298, USA
 May be reproduced by subscribing institution for internal distribution only.