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ANEXO

www.ecri.org . Printed from Health Devices Alerts on Tuesday, May 29, 2018 Page 1

[High Priority] - A30330 : Roche—cobas e 801 Module Plungers: Syringe Function May Be Impaired

Medical Device Ongoing Action

Published: Thursday, May 24, 2018

UMDNS Terms:

- Analyzers, Laboratory, Clinical Chemistry/Immunoassay [20821]

Product Identifier:
[Capital Equipment]

Product	Hitachi High-Technologies Corp (Manufacturer) Model	Part No.	Serial No.	Item No.
Syringe Assembly Plungers	cobas e 801 Module PreWash (PW), cobas e 801 Module R1, cobas e 801 Module R2	798-3203	None listed	None listed
Modules	cobas e 801	None listed	1601-01 to 18E6-10	07682913001

Geographic Regions: Argentina, Australia, Brazil, Canada, Chile, Colombia, Europe, Hong Kong, Japan, Korea, Malaysia, Middle East, New Zealand, South Africa, Thailand, Turkey, U.K., U.S., Vietnam

Manufacturer(s): Hitachi High-Technologies Corp (Manufacturer)1-24-14 Nishi-Shinbashi, 105-8717 Minato-ku Tokyo, Japan
Roche Diagnostics GmbH (Authorized Representative)Sandhoferstrasse 116, Mannheim, D-68305, Germany

Suggested Distribution: Clinical/Biomedical Engineering, Clinical Laboratory/Pathology

Problem: In a March 2018 Urgent Field Safety Notice letter posted by the U.K. Medicines and Healthcare Products Regulatory Agency (MHRA) and the German Federal Institute for Drugs and Medical Devices (BfArM), Roche states that it has received four complaints of a problem with the syringe plunger on a cobas e 801 module PreWash (PW) syringe assembly. The problem may also occur on the cobas e 801 module R1 or R2 reagent syringe assembly. Depending on the affected syringe, R1, R2, or prewash liquids will not be pipetted correctly, as follows:

- (1) If the PreWash syringe plunger is broken, the cobas e 801 module cannot detect it and no alarm is issued. In this instance, the prewash step required by most assays will not be performed and an effect on measurement results cannot be excluded. The "Important Information" section in the [letter](#) outlines assays which are not affected by the problem. All other assays run on the cobas e 801 module may be affected.
- (2) If the R1 or R2 syringe plunger is broken, a system alarm "345-1 Warning Abnormal Low Signal" is generated and data alarm (-SigL) is flagged; consequently, no measurement results will be generated.

Depending on how the seal pieces are assembled, the plunger may be mounted in a tilted position, which applies mechanical stress to the plunger. In the worst case scenario during the cobas e 801 initialization steps, when the largest strokes happen, the plunger may crack leading to an impaired syringe function.

Action Needed:

Identify any affected product in your inventory. If you have affected product, verify that you have received the March 2018 Urgent Field Safety Notice letter from Roche. As a final solution, a mandatory installation of a syringe plunger guide for the R1, R2, and PW syringes is required for cobas e 801 modules in the field, from serial numbers 1601-01 to 18E6-10. The hardware modification kit is already available in the warehouse. The installation of the kit requires a few minutes, if the instrument is available in [Standby] mode. The kit ensures the correct positioning and movement of the above syringe plungers and thus prevents the problem. Your Roche Diagnostics local representative will contact your facility to schedule the installation of the modification kit. The installation process is ongoing. If the described problem occurs with an R1 or R2 syringe plunger:

- (1) A system alarm "345-1 Warning Abnormal Low Signal" is generated.
- (2) A data alarm (-SigL) is flagged.
- (3) No measurement results are generated.
- (4) Contact Roche's customer support department.

If the described problem occurs with the PW syringe plunger:

- (1) Refer to Table 1 in the [letter](#) for a list of cobas e 801 assays that are not affected because they do not require the prewash step.
- (2) Refer to Table 2 in the [letter](#) for a list of cobas e 801 assays with controls within range and samples within control range. These assays require prewash steps. The effect on measurement results by omitting prewash steps is assay-dependent. Studies were performed to investigate the impact of a missing prewash step on controls and samples.
- (3) Refer to Table 3 in the [letter](#) for a list of cobas e 801 prewash assays. These assays require prewash steps. The effect on measurement results by omitting prewash steps is assay-dependent. For these assays, the effect on measurement results cannot be excluded.

If it is not possible to have the above mentioned final solution (modification kit) in a timely manner, follow temporarily one of the below listed options (workarounds):

Option 1 (workaround):

This workaround (option 1) can be applied by the operators without support from a Roche field service representative. Measure only assays that do not require prewash steps and assays which are not affected by the missing prewash step.

The lists are found under "Important Information", Table 1 and Table 2 in the [letter](#).

The assays listed in Table 3 under "Important Information" must NOT be measured until installation of the modification kit has been completed. These assays have to be masked. The procedure for test masking is available in the cobas 8000 Operator's Manual, under Operation > Order and results > Processing samples > Masking tests.

Option 2 (workaround):

This workaround (option 2) requires intervention from your Roche local field service representative for the necessary configuration of the data manager. Consequences on the laboratory workflow will need to be considered.

Visually inspect the functionality of the PreWash syringe plunger before manually releasing results to Host/LIS, until the final solution is applied.

The automatic data transfer function of patient sample results on data manager needs to be disabled, so patient results measured on cobas e 801 module (prewash assays only) will not be automatically transferred to Host/LIS. Instead, these measurement results must be manually released to Host/LIS after visually checking the functionality of the PW syringe plunger.

Option 2 - Procedure:

- (1) Define the time interval to visually inspect the PreWash syringe plunger and manually validate results (for prewash assays) on the data manager to be sent to Host/LIS.

©2018 ECRI Institute
5200 Butler Pike, Plymouth Meeting, PA 19462-1298, USA
May be reproduced by subscribing institution for internal distribution only.

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www.ecri.org . Printed from *Health Devices Alerts* on Tuesday, May 29, 2018 Page 2

(2) Start sample processing, and after the abovedefined time interval, stop the operation with the [5.Stop] button and waituntil the system goes to [Standby] mode.

(3) Visually check the PreWash syringe plunger movement, according to the attachment:

"Att. 01_FSN_PreWash Plunger check during System Prime"

(4) If the PreWash syringe plunger is visually working as expected, manually validate the transfer of patient results measured on cobas e 801 from data manager to Host/LIS. During the manual validation of patient results, it is recommended that the system is kept in [Standby] mode and no new samples are loaded.

(5) Repeat Step 2 through 4 until the end of routine operation.

If a damaged PreWash syringe plunger is identified (see the above procedure, step 3), stop using the cobas e 801 module and immediately contact Roche's customer support department. In addition, do not validate the measured results and repeat the measurement of concerned samples.

Notify all relevant personnel at your facility of the information in the letter, and forward a copy of the letter to any facility to which you have further distributed affected product.

For Further Information:

Roche

Website: [Click here](#)

References:

- Great Britain. Medicines and Healthcare Products Regulatory Agency. Hitachi: cobas 8000 modular analyzer cobas e801 module [online]. London: Department of Health; 2018 Apr 3 [cited 2018 May 22]. (Field safety notice; reference no. 2018/003/023/478/028). Available from Internet: [Click here](#) .
- Germany. Federal Institute for Drugs and Medical Devices. Urgent Field Safety Notice for cobas 8000 by Hitachi High-Technologies Corporation [online]. 2018 Apr 16 [cited 2018 May 22]. Available from Internet: [Click here](#) .

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#) .

Source(s):

- 2018 May 22. MHRA FSN. 2018/003/023/478/028 [Download](#)
- 2018 May 22. MHRA FSN. SBN-CPS-2018-002 [Download](#)
- 2018 May 22. BfArM (Germany). 03589/18 [Download](#)
- 2018 May 22. BfArM (Germany). SBN-CPS-2018-002 [Download](#)
- 2018 May 22. Manufacturer. Manufacturer confirmed information contained in source material