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ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Monday, April 3, 2017 Page 1

[High Priority] - A28218 : GE— Vivid Ultrasound Systems: User May Select Incorrect Patient if DICOM Worklist Search Response Time Is Slow
Medical Device Ongoing Action

Published: Thursday, March 16, 2017
 Last Updated: Monday, April 3, 2017

UMDNS Terms:

- Scanning Systems, Ultrasonic, General-Purpose [15976]
- Scanning Systems, Ultrasonic, Cardiac [17422]
- Scanning Systems, Ultrasonic, Vascular [15957]

Product Identifier:

Vivid Ultrasound Systems:	Software Versions:
E80, E90, E95	201, revision 54.0, 61.0
S60, S70	201, revision 55.0, 63.0
S60N, S70N	201, revision 63.0

[Capital Equipment]

Geographic Regions: Worldwide

Manufacturer(s): GE Healthcare9900 Innovation Dr, Wauwatosa, WI 53226, United States

Suggested Distribution: Cardiology/Cardiac Catheterization Laboratory, Clinical/Biomedical Engineering, Diagnostic Imaging, Information Technology

Problem:

In a March 10, 2017, Urgent Medical Device Correction letter submitted by an ECRI Institute member hospital, GE states that the user may incorrectly select a patient in the above systems when the DICOM Worklist search response time is slow, potentially leading to misdiagnosis if the user does not recognize that images from the intended patient are stored under the incorrect patient after the exam. GE also states that it has received no reports of injuries as a result of this problem.

Action Needed:

Identify any affected systems in your inventory. If you have affected systems, verify that you have received the March 10, 2017, Urgent Medical Device Correction letter from GE. To reduce the likelihood of this problem occurring, reduce the size of the DICOM Worklist search by setting the "max results" to 75 in the Worklist dataflow config dialog and/or configuring the DICOM Worklist server to automatically delete "Scheduled Procedure Steps" that have been performed on the ultrasound scanner. If the latter option is not available on your DICOM Worklist server, manually clear all "Scheduled Procedure Steps" that have been performed from the DICOM Worklist on a regular basis. Consult your facility's IT department for assistance with this option. After selecting a patient in the DICOM Modality Worklist, always review and confirm that the patient demographics shown on the ultrasound scanner screen during the exam are from the intended patient. If the incorrect and/or intended patient already exist in the patient archive, a Patient Match dialog will display to inform the user about a patient information inconsistency. A GE representative will contact your facility to arrange for a correction at no cost.

For Further Information:

GE
 Tel: (800) 437-1171
 Website: [Click here](#)

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

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