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ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Wednesday, April 12, 2017 Page 1

[High Priority] - A28324 : SCHILLER—FRED easy Defibrillator Batteries: May Fail When Inserted into Defibrillator

Medical Device Ongoing Action

Published: Monday, April 3, 2017

UMDNS Terms:

- Defibrillators, External, Automated [17116]
- Batteries [16640]

Product Identifier:

Batteries for use with FRED easy Automated External Defibrillators (AEDs) [Consumable, Capital Equipment]
Battery Reference No. 4-07-0001; Battery Lot Nos.: HYB 16-0056, HYB 16-0057, HYB 16-0058, HYB 16-0059
Batteries sold between October 1, 2016, and January 15, 2017

Geographic Regions: Asia, Europe, Latin America

Manufacturer(s): SCHILLER AG Altgasse 68, CH-6341, Baar, Switzerland

Suggested Distribution: Cardiology/Cardiac Catheterization Laboratory, Clinical/Biomedical Engineering, Critical Care, Emergency/Outpatient Services, Nursing, Home Care, EMS/Transport, Materials Management

Problem:

In a January 17, 2017, Field Safety Notice letter posted by the German Federal Institute for Drugs and Medical Devices (BfArM), SCHILLER states that the above batteries may cause a "BATTERY EMPTY" failure when inserted into the above AEDs. SCHILLER also states that batteries with this fault cannot be used and that the fault is detected automatically by the defibrillator when the battery is inserted into the device.

Action Needed:

Identify any affected product in your inventory. For an image of the location of the reference number and lot number on the battery packaging, see the image in the Affected Devices section of the [letter](#). If you have affected product, verify that you have received the January 17, 2017, Field Safety Notice letter from SCHILLER. Test affected batteries by inserting them into your AED (for instructions on how to do this, see section 3.1 "Inserting the battery" in the FRED easy user manual). The unit will perform its self-test. At the end of the self-test, the unit will shut down. If the green LED flashes, the battery is free from defects and can be used without limitation. For an image of the green LED, see the [letter](#). If, after the self-test, the green LED remains off and the display reads "EMPTY BATTERY>REPLACE BATTERY" (for an image of this display, see the [letter](#)), the battery is defective. Contact SCHILLER to arrange for defective battery replacement.

For Further Information:

SCHILLER

Website: [Click here](#)

References:

- Germany. Federal Institute for Drugs and Medical Devices. Urgent field safety notice for defibrillator battery FRED easy, Schiller Medical SAS [online]. 2017 Mar 28 [cited 2017 Mar 29]. Available from Internet: [Click here](#).

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2017 Mar 29. BfArM (Germany). 01529/17 [Download](#)
- 2017 Mar 29. BfArM (Germany). [Download](#)
- 2017 Apr 3. Manufacturer. SCHILLER confirmed the information provided in the source material.