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ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Wednesday, April 12, 2017 Page 1

[High Priority] - A28314 : Draeger—Infinity Monitors Used with Scio Modules: May Not Generate Certain Alarms

Medical Device Ongoing Action

Published: Monday, April 3, 2017

UMDNS Terms:

- Monitors, Telemetric, Physiologic [13987]

Product Identifier:

Infinity Monitors: (1) Delta, (2) Delta XL, (3) Kappa used with the following Modules: (4) Scio, (5) Scio Four, (6) Scio Four Oxi Plus, (7) Scio Four Plus [Capital Equipment]
Software Version VF10.0

Geographic Regions: Worldwide

Manufacturer(s): Draeger Medical Deutschland GmbH Moislinger Allee 53-55, Luebeck, D-23558, Germany

Suggested Distribution: Clinical/Biomedical Engineering, Critical Care, Emergency/Outpatient Services, Nursing, OR/Surgery, Pulmonology/Respiratory Therapy, Information Technology, NICU

Problem: In a March 2017 Important Safety Notice letter posted by the U.K. Medicines and Healthcare Products Regulatory Agency (MHRA), Draeger states that the above monitors with software VF10.0 when combined with the referenced Scio modules will not generate the following alarms even though the level measured by the Scio module is outside of the set alarm limit: (# represents the set alarm limit): $iO_2 > \#$, $iO_2 < \#$, $i[agent] > \#$, $i[agent] < \#$, Inspired MAC High, $etN_2O > 82\%$. Draeger states that although the alarms are not generated, the parameter readings are not affected and will be displayed properly.

Action Needed:

Identify any affected product in your inventory. If you have any affected product, verify that you have received the March 2017 Important Safety Notice/Medical Device Recall letter and Reply and Order Card from Draeger. If you are using the above referenced Scio modules in your facility, do not combine them with Delta, Delta XL or Kappa monitors running software VF10.0 unless you can ensure adequate patient surveillance. A downgrade to software version VF9.1 is required as a short-term measure. If you have the above modules, contact Draeger using the attached Reply and Order card. Your Draeger sales and service representative will contact your facility to install the downgrade to VF9.1 at no charge. Inform all relevant personnel of the information in the Urgent Field Safety Notice letter. If you do not have any modules, no immediate action is required. Draeger is developing a new software version VF10.1. Once the VF10.1 Software is released, your Draeger sales and service representative will contact your facility to arrange to install the update to VF10.1 at no charge. Delta Family monitors that were downgraded because of this problem will also be upgraded to VF10.1 software.

For Further Information:

Draeger
Website: [Click here](#)

References:

- Great Britain. Medicines and Healthcare Products Regulatory Agency. Draeger: Delta/Delta XL/Kappa [online]. London: Department of Health; 2017 Mar 27 [cited 2017 Mar 29]. (Field safety notice; reference no. 2017/003/021/291/007). Available from Internet: [Click here](#).

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2017 Mar 29. MHRA FSN. 2017/003/021/291/007 [Download](#)
- 2017 Mar 29. MHRA FSN. [Download](#)
- 2017 Mar 30. Manufacturer. Manufacturer confirmed source documents