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### ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Wednesday, March 22, 2017 Page 1

#### [Normal Priority ] - A28134 : Tosoh—Model HLC-723G8 Glycohemoglobin Analyzers: May Yield Negatively Biased HbA1c Results Medical Device Ongoing Action

Published: Wednesday, March 1, 2017  
Last Updated: Thursday, March 2, 2017

#### UMDNS Terms:

- IVD Test Reagent/Kits, Clinical Chemistry, Protein, Glycated Hemoglobin [18904]
- Analyzers, Laboratory, Blood, Glycated Hemoglobin [17109]

**Product Identifier:** Model HLC-723G8 Glycohemoglobin Analyzers [Capital Equipment]  
All active analyzers  
682 units distributed

**Geographic Regions:** Asia, Brazil, Canada, Chile, Colombia, Costa Rica, Dominican Republic, Ecuador, Europe, Guatemala, Nicaragua, Panama, Peru, Puerto Rico, Uruguay, U.S.

**Manufacturer(s):** Tosoh Hi-Tec 1-37 Fukugawa minami-machi, Syunan-shi, Yamaguchi, 746-0042, Japan

**Suggested Distribution:** Clinical/Biomedical Engineering, Clinical Laboratory/Pathology

#### Summary:

On March 1, 2017, the manufacturer confirmed the information in the source material and verified additional geographic regions of distribution, bolded above.

**Problem:** FDA's Center for Devices and Radiological Health (CDRH) states that the above analyzers may yield HbA1c values with a negative bias for patients with hemoglobin variants C, D, E, and S with higher levels of HbA1c (e.g., 9%), but not at the diagnostic medical decision level of 6.5%. FDA's CDRH also states that the manufacturer initiated a recall by Recall Notification letter on November 7, 2016.

#### Action Needed:

Identify any affected analyzers in your inventory. If you have affected analyzers, verify that you have received the Recall Notification letter from Tosoh. FDA's CDRH states that users should review the chromatography to determine whether a variant peak was detected and, if so, retest the sample using an alternative method.

#### For Further Information:

Susan H. Koss, Tosoh  
Tel.: (614) 317-1909  
Website: [Click here](#)

#### References:

- United States. Food and Drug Administration. Center for Devices and Radiological Health. Class 2 device recall Tosoh Bioscience HLC723GB analyzer [online]. 2017 Feb 17 [cited 2017 Feb 27]. Available from Internet: [Click here](#).

#### Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

#### Source(s):

- 2017 Feb 27. FDA CDRH Database. Class II Z-1225-2017 [Download](#)
- 2017 Mar 1. Manufacturer. The manufacturer confirmed the information in the source material.