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### ANEXO 1

www.ecri.org . Printed from *Health Devices Alerts* on Thursday, March 23, 2017 Page 1

#### [High Priority ] - A28125 : ELLA-CS—Stent Delivery System Fixation-Gastric Balloons Contained in Danis Procedure Packs: Safety Valve May Exhibit Unintended Movement Medical Device Ongoing Action

Published: Monday, March 6, 2017  
Last Updated: Monday, March 20, 2017

#### UMDNS Terms:

- Procedure Trays, Surgical [17168]
- Stents, Esophageal, Metallic, Self-Expanding [20444]

**Product Identifier:** Stent Delivery System Fixation-Gastric Balloons contained in Danis Procedure Packs [Consumable]  
Procedure Pack Reference Nos.: 019-08S-25-135, 019-08S-25-135-B

**Geographic Regions:** Australia, Austria, Bahrain, Belarus, Belgium, Bolivia, Brazil, Canada, Chile, Colombia, Czech Republic, Denmark, Ecuador, Egypt, Finland, France, Germany, Georgia, Hong Kong, India, Ireland, Italy, Israel, Kuwait, Malaysia, Moldova, Myanmar, The Netherlands, New Zealand, Norway, Panama, Peru, Poland, Russian Federation, Saudi Arabia, Slovakia, Slovenia, South Africa, Spain, Sweden, Switzerland, Turkey, U.K., Vietnam

**Distributor(s):** • UK Medical Ltd A CareFusion Company Albreda House, Lydgate Lane, Sheffield S105FH England (U.K. distributor)

**Manufacturer(s):** ELLA-CS s r oMilady Horakove 504/45, Trebes 400 Hradec Kralove, Czech Republic

**Suggested Distribution:** Emergency/Outpatient Services, OR/Surgery, Gastroenterology, Materials Management

#### Problem:

In a February 14, 2017, Urgent Field Safety Notice letter posted by the U.K. Medicines and Healthcare Products Regulatory Agency (MHRA) and the German Federal Institute for Drugs and Medical Devices (BfArM), ELLA-CS states that the safety valve of the fixation/gastric balloons of the above Danis stent delivery systems may exhibit unintended movement, potentially resulting in the following:

- (1) Inability of the physician to correctly place the delivery system because it is impossible to insufflate the fixation-gastric balloon.
- (2) The physician may need to use another Danis stent or other method to stop acute esophageal variceal bleeding, potentially prolonging the bleeding and/or the procedure.

#### Action Needed:

Identify any affected product in your inventory. If you have affected product, verify that you have received the February 14, 2017, Urgent Field Safety Notice letter and Acknowledgment of Receipt Form from ELLA-CS. The firm recommends the following actions:

- (1) Check the position of the safety valve (refer to the figure in the [letter](#)).
- (2) If the safety valve does not cover the opening designed for spontaneous deflation of the fixation-gastric balloon, return it to the appropriate position to allow performance of the stent implantation with the fixation of the delivery system in the cardia.
- (3) If the safety valve does not cover the opening designed for spontaneous deflation of the fixation-gastric balloon and is not returned to the appropriate position, the Danis stent may be introduced using conventional techniques for implantation of esophageal stents (i.e., under fluoroscopic or combined endoscopic and fluoroscopic guidance).

Inform all relevant personnel at your facility of the information in the Urgent Field Safety Notice letter. Complete the Acknowledgment of Receipt Form, and return it to ELLA-CS using the instructions on the form.

#### For Further Information:

ELLA-CS  
E-mail: [vigilance@ellacs.eu](mailto:vigilance@ellacs.eu)  
Website: [Click here](#)

#### References:

- Great Britain. Medicines and Healthcare Products Regulatory Agency. ELLA-CS: SX-ELLA stent Danis procedure pack (basic) [online]. London: Department of Health; 2017 Feb 20 [cited 2017 Feb 28]. (Field safety notice; reference no. 2017/002/015/291/004). Available from Internet: [Click here](#).
- Germany. Federal Institute for Drugs and Medical Devices. Urgent field safety notice for Danis procedure pack/Basic, ELLA-CS s.r.o. [online]. 2017 Mar 2 [cited 2017 Mar 2]. Available from Internet: [Click here](#).

#### Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

#### Source(s):

- 2017 Feb 28. MHRA FSN. 2017/002/015/291/004 [Download](#)
- 2017 Feb 28. MHRA FSN. ELLA-CS Reference No. FSN-ELLA-2017-001 [Download](#)
- 2017 Mar 2. BfArM (Germany). 01508/17 [Download](#)
- 2017 Mar 2. BfArM (Germany). ELLA-CS Reference No. FSN-ELLA-2017-001 [Download](#)
- 2017 Mar 6. Manufacturer. The manufacturer confirmed the information in the source material.