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High Priority

Published: Thursday, January 11, 2018

Beckman Coulter—Access 2 and UniCel DxC 600i SYNCHRON Access Clinical Systems: Skin Burn May Occur if Tool Contacts Reagent Storage Interface Board When Pipettor Alignment Is Performed

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Geographic Regions: (Impact in additional regions has not been identified or ruled out at the time of this posting), Argentina, Canada, Chile, China, Colombia, Costa Rica, France, Georgia, Germany, Ghana, Greece, Hong Kong, India, Italy, Japan, Kazakhstan, Mexico, Morocco, Myanmar, Nepal, Philippines, Poland, Russia, Rwanda, Saudi Arabia, Slovakia, South Africa, South Korea, Thailand, Turkey, U.K., U.S., Vietnam

Manufacturer(s): Beckman Coulter Inc [362269], 250 S Kraemer Blvd, Brea, CA 92821-6232, United States

Problem:

FDA's Center for Devices and Radiological Health (CDRH) states that minor skin burns may occur if the tool contacts the reagent storage interface board (B78817) when a pipettor alignment is performed on the above systems. FDA's CDRH also states that the manufacturer initiated a recall by Urgent Medical Device Recall letter dated November 6, 2017. The manufacturer has not confirmed the information provided in the source material.

Action Needed:

Identify any affected systems in your inventory. If you have affected systems, verify that you have received the November 6, 2017, Urgent Medical Device Recall letter from Beckman Coulter. Beckman Coulter states that you may continue to operate affected systems normally. If a device on the instrument needs alignment, contact the Beckman Coulter technical support department. The alignment tool is to be used only by Beckman Coulter service personnel and Beckman Coulter trained biomedical engineers. A Beckman Coulter service representative will contact your facility to schedule a service visit. Notify all relevant personnel at your facility of the information in the letter, and retain a copy of the letter with your records. Respond to Beckman Coulter using the instructions in the letter.

For Further Information:

Beckman Coulter
 Website: [Click here](#)

References:

United States:

- Food and Drug Administration. Center for Devices and Radiological Health. Class 3 device recall Access 2 immunoassay system [online]. 2018 Jan 3 [cited 2018 Jan 10]. Available from Internet: [Click here](#).
- Food and Drug Administration. Center for Devices and Radiological Health. Class 3 device recall Access 2 immunoassay system [online]. 2018 Jan 3 [cited 2018 Jan 10]. Available from Internet: [Click here](#).
- Food and Drug Administration. Center for Devices and Radiological Health. Class 3 device recall Access 2 immunoassay system [online]. 2018 Jan 3 [cited 2018 Jan 10]. Available from Internet: [Click here](#).

Comment:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Suggested Distribution: Clinical Laboratory/Pathology, Clinical/Biomedical Engineering