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ANEXO

www.ecri.org . Printed from *Health Devices Alerts* on Wednesday, June 13, 2018 Page 1

[High Priority] - A30604 01 : GE/Silex: Various MAC ECG Systems and Wireless Device Servers: May Be Vulnerable to Improper Authentication and OS Command Injection [Update]
 Medical Device Ongoing Action

Published: Monday, June 4, 2018

UMDNS Terms:

- Information Systems, Data Management, Cardiology, Electrocardiography [22499]
- Network Servers [21973]

Product Identifier:
 [Capital Equipment]

Product	GE Healthcare Model	
Resting ECG Analysis Systems	MAC 3500, MAC 5000, MAC 5500, MAC 5500 HD	
Wireless ECG Communication Modules	MobileLink	
Product	Silex Technology Inc Model	Version
Wireless Device Servers	GEH-500	<= 1.54
	SX-500	All
	GEH-SD-320AN	<= GEH-1.1
	SD-320AN	<= 2.01

Geographic Regions: (Impact in specific regions has not been identified or ruled out at the time of this posting), Worldwide

Manufacturer(s): GE Healthcare 9900 Innovation Dr, Wauwatosa, WI 53226, United States
 Silex Technology Inc 167 West 7065 South, Midvale, UT 84047, United States

Suggested Distribution: Cardiology/Cardiac Catheterization Laboratory, Clinical/Biomedical Engineering, Critical Care, Emergency/Outpatient Services, Information Technology, EMS/Transport

Summary:

Update Reason: GE and Silex provide updated firmware image for the GEH-SD-320AN system. This Alert provides additional information based on a May 31, 2018, Advisory revision regarding [Alert A30604](#). The updated firmware image is available for download [here](#).

Problem:

In a May 8, 2018, Advisory, the U.S. Industrial Control Systems Cyber Emergency Response Team (ICS-CERT) states that the above systems may be susceptible to the following vulnerabilities:

- Improper authentication (CWE-287)
 - Authentication is not verified when making certain POST requests, which may allow attackers to modify system settings.
- Improper neutralization of special elements used in an OS command (OS Command Injection) (CWE-78)
 - A system call parameter is not properly sanitized, which may allow remote code execution.

The manufacturers have not confirmed the information provided in the source material.

Action Needed:

The following actions are those listed in [Alert A30604](#). Identify any affected systems in your inventory. If you have affected systems, verify that you have reviewed the May 8, 2018, [ICS-CERT Advisory](#). GE and Silex recommend the following mitigations:

- For MobileLink and SX-500 systems, enable the "update" account within the web interface, which is not enabled by default. Set the secondary password for the "update" account to prevent unauthenticated changes to the device configuration.
- For MobileLink and GEH-SD-320AN systems, Silex Technology and GE Healthcare have produced an updated firmware image for the GEH-SD-320AN system, which will be made available for download from GE upon completion of testing.

GE will post information pertaining to enabling the "update" account and download of new firmware [here](#). The firmware update for SD-320AN is separate from GEH-SD-320AN and will be available for download from Silex Technology at a future date. This update does not pertain to the listed GEH device. Contact Silex Technology for more information regarding download and application of this new firmware. National Cybersecurity and Communications Integration Center (NCCIC) recommends the following actions to minimize the risk of exploitation of these vulnerabilities:

- Minimize network exposure for all control system devices and/or systems, and ensure that they are not accessible from the Internet.
- Locate control system networks and remote devices behind firewalls, and isolate them from the business network.
- When remote access is required, use secure methods, such as Virtual Private Networks (VPNs), recognizing that VPNs may have vulnerabilities and should be updated to the most current version available. Also recognize that a VPN is only as secure as the connected

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www.ecri.org . Printed from *Health Devices Alerts* on Wednesday, June 13, 2018 Page 2

devices.

Perform proper impact analysis and risk assessment before deploying these defensive measures. Additional mitigation guidance and recommended practices are publicly available in the NCCIC Technical Information Paper, [ICS-TIP-12-146-01B--Targeted Cyber Intrusion Detection and Mitigation Strategies](#), which is available for download from the [ICS-CERT website](#). NCCIC also provides a section for [control systems security recommended practices](#). Several recommended practices are available for reading and download, including [Improving Industrial Control Systems Cybersecurity with Defense-in-Depth Strategies](#).

For Further Information:

NCCIC

E-mail: NCCICCUSTOMERSERVICE@hq.dhs.gov

Tel.: (888) 282-0870

Website: [Click here](#)

GE

Website: [Click here](#)

Silex

Website: [Click here](#)

References:

United States:

- Department of Homeland Security. Industrial Control Systems Cyber Emergency Response Team. Advisory. Silex Technology SX-500/SD-320AN or GE Healthcare MobileLink [online]. 2018 May 8 [cited 2018 May 14]. Available from Internet: [Click here](#).
- Department of Homeland Security. Industrial Control Systems Cyber Emergency Response Team. Advisory. Silex Technology SX-500/SD-320AN or GE Healthcare MobileLink (Update A) [online]. 2018 May 31 [cited 2018 Jun 1]. Available from Internet: [Click here](#).

Comments:

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our [HDA Format Guide](#).

Source(s):

- 2018 Jun 1. ICS-CERT Advisory: ICSMA-18-128-01 [Download](#)