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### ANEXO

www.ecri.org . Printed from *Health Devices Alerts* on Friday, June 29, 2018 Page 1

**[High Priority ] - A30849 : Natus—Xitek NeuroWorks Software: May Be Susceptible to Cybersecurity Vulnerabilities  
Medical Device Ongoing Action**

**Published:** Tuesday, June 19, 2018

**UMDNS Terms:**

- Software, Physiologic Monitoring, Electroencephalography [26714]

**Product Identifier:**  
[Capital Equipment]

Product	Natus Medical Inc Model	Software Version
Electroencephalography (EEG) System Software	Xitek NeuroWorks	8

**Geographic Regions:** (Impact in specific regions has not been identified or ruled out at the time of this posting), Worldwide

**Manufacturer(s):** Natus Medical Inc 6701 Koll Center Pkwy, Pleasanton, CA 94566, United States

**Suggested Distribution:** Clinical/Biomedical Engineering, Critical Care, Information Technology, Neurology, Sleep Laboratory

**Problem:**

In a June 14, 2018, Advisory, the U.S. Industrial Control Systems Cyber Emergency Response Team (ICS-CERT) states that the above software may be vulnerable to Stack-Based Buffer Overflow and Out-of-Bounds Read. ICS-CERT also states that successful exploitation of these vulnerabilities require access to the Natus customer network, and could crash the device being accessed; a buffer overflow condition may allow remote code execution. The potential vulnerabilities are as follows:

- Out-of-Bounds Read (CWE-125)
  - A specially-crafted packet may be able to cause an out-of-bounds read, which may result in a denial-of-service condition.
- Stack-Based Buffer Overflow (CWE-121)
  - An attacker may cause a buffer overflow by sending a specially-crafted packet to the affected product while the product attempts to open a file requested by the client.
  - A specially-crafted packet received during the execution of certain commands can cause memory to be overwritten in a way that could allow an attacker to take control of the program.
  - An error in the way the program parses data structures may allow an attacker to take control of the system by sending it a specially-crafted packet.
  - A specially-crafted packet takes advantage of the way the program parses data structures and may cause a buffer overflow, which may allow remote execution of arbitrary code.

The manufacturer has not confirmed the information provided in the source material.

**Action Needed:**

Identify any affected software in your inventory. If you have affected software, verify that you have reviewed the June 14, 2018, [ICS-CERT Advisory](#). Natus recommends installation of NeuroWorks/SleepWorks 8.5 GMA 3, a software update with security enhancements to address the vulnerabilities identified in NeuroWorks/SleepWorks 8. A free software update to NeuroWorks/SleepWorks 8.5 GMA 3 is available to users using NeuroWorks/SleepWorks Version 8.0, 8.1, 8.4, or 8.5. Natus recommends installing this update on affected systems as quickly as possible. The National Cybersecurity and Communications Integration Center (NCCIC) recommends that users take defensive measures to minimize the risk of exploitation of the above vulnerabilities. Users should:

- Minimize network exposure for all control system devices and/or systems, and ensure that they are not accessible from the Internet.
- Locate control system networks and remote devices behind firewalls, and isolate them from the business network.
- When remote access is required, use secure methods, such as virtual private networks (VPNs), recognizing that VPNs may have vulnerabilities and should be updated to the most current version available. Recognize that VPN is only as secure as the connected devices.

Additional mitigation guidance and recommended practices are publicly available on the [ICS-CERT website](#) in the Technical Information Paper, [ICS-TIP-12-146-01B--Targeted Cyber Intrusion Detection and Mitigation Strategies](#).

**For Further Information:**

NCCIC  
 Tel.: (888) 282-0870  
 E-mail: [NCCICCUSTOMERSERVICE@hq.dhs.gov](mailto:NCCICCUSTOMERSERVICE@hq.dhs.gov)  
 Website: [Click here](#)  
 Natus  
 Website: [Click here](#)

**References:**

- United States. Department of Homeland Security. Industrial Control Systems Cyber Emergency Response Team. Advisory. Natus Xitek NeuroWorks [online]. 2018 Jun 14 [cited 2018 Jun 19]. Available from Internet: [Click here](#).

**Comments:**

©2018 ECRI Institute  
 5200 Butler Pike, Plymouth Meeting, PA 19462-1298, USA  
 May be reproduced by subscribing institution for internal distribution only.

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www.ecri.org . Printed from *Health Devices Alerts* on Friday, June 29, 2018 Page 2

- This alert is a living document and may be updated when ECRI Institute receives additional information. In circumstances in which we determine that it is appropriate for customers to repeat their review of an issue (e.g., when additional affected product has been identified), we will post a separate update alert. In other cases, we may add information, such as additional commentary, recommendations, and/or source documents, to the original alert. For additional information regarding the format of this alert, refer to our HDA Format Guide.

**Source(s):**

- 2018 Jun 19. Member Hospital. ICS-CERT Advisory [Download](#)